

Automatic Blood Pressure Monitor, Wrist Model

HCPCS: A4670-1,A4670



FIGURE A

- Large easy-to-read digital display
- 120 memory recall
- Date & time display
- Pulse Rate Indicator
- Auto-off function
- Lightweight & easy to use
- Quiet inflation method
- Wrist model

PRODUCTS	Description
BP2116	BP Monitor Wrist Auto 1 EA
SPECIFICATIONS	Description
Auto Power Off	3 Mnutes After Last Reading
Measurement Method	Oscillometric
Measurement Range	Pressure: 0~300mmHg +/-3mmHg
Memory Capacity	120 Memories with Date & Time
Operation Environment	50F tp 104F, <85%RH
Power Supply	2 x AAA Batteries
Warranty	2 Years