


## Drive Regulators

HCPCS:



- Drive regulators feature a protective contents gauge and a built-in pressure relief valve for safety
- Mini regulators offer the same function as standard regulators while facilitating an easier fit into oxygen carry bags
- Easy-to-operate single selector knob
- Lightweight and durable design
- Meet or exceed accuracy standards for ASTM, American National and CGA
- Inlet Pressure: 300 - 3,000 PSI
- Outlet pressure: 50 PSI
- UL Approved

PRODUCTS	Description
18300GN	Drive Regulators, EA
18301G	Drive Regulators, BX
18301GM	Drive Regulators, BX
18302GDEL	Drive Regulators, EA
18302GMN	Drive Regulators, EA
18302GN	Drive Regulators, EA
18303G	Drive Regulators, BX
18303GM	Drive Regulators, BX
18304GMN	Drive Regulators, EA
18304GN	Drive Regulators, EA
18307G	Drive Regulators, EA
SPECIFICATIONS	Description
Warranty	Limited 3 Year
ACCESSORIES	Description
SEAL3000	Seal Washer, CHAD O2 Cnsvrs&Regs,10/pkg, PKG
OM-900-FEC	E-Clip for Conserver Handles, 10/bag, BAG