

Roscoe Mobility

E-Series Rollator

The Roscoe E-Series Rollator will meet all your patient needs with a focus on style and affordability. Easy-to-operate ergonomic hand brakes, zippered nylon pouch complete with a flip-up curved padded backrest add to the comfort of the Roscoe E-Series Rollator. Folds up flat for easy transport and storage. The Roscoe E-Series is designed to meet all your fundamental rollator needs and will meet a variety of patient profiles.



Burgundy
Finish



Blue
Finish



Product Features


- Lightweight aluminum frame
- 6" wheels with wheel locks
- Easy-to-operate, ergonomic loops-style brakes
- Padded seat with easy-folding mechanism
- Flip-up Curved padded backrest
- Zippered compartment under seat for storage of personal items
- One year warranty
90 days on all other parts

Product Specifications

- Seat Measurements: 13 3/4" (W) x 11 1/2" (D)
- Seat-to-Floor Height: 21"
- Adjustable Handle Height: 31" - 36"
- Unit Width: 22 3/4"
- Product Weight: 13 lb.
- Weight Capacity: 300 lb.
- HCPCS Code: E0143*
- HCPCS Code: E0156 for seat



Locking Loop Style Brakes



Flip-Up Curved Padded Backrest



Flip-up Seat with Zippered Pouch

Item #	Description	Color	Qty
30164	E-Series Rollator with Padded Seat, 6-inch Wheels	Burgundy	each
30166	E-Series Rollator with Padded Seat, 6-inch Wheels	Blue	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com