

ENDURANCE[®] HEIGHT ADJUSTABLE HAND BED RAILS

Give your customers the “Helping Hand” they need to get in and out of bed safely



P1410 HAND BED RAIL

- 20" wide handle has ¼" thick foam tubing to provide a secure grip
- Slides up to 43" under mattress with easy adjustments and has non-slip sleeves to keep it in place with adjustable strap to secure to frame
- Breaks down for travel
- Made of chrome plated steel

SPECIFICATIONS

User Weight Limit: 300 lbs.
 Height Adj: 17" - 21"



Introducing the **NEW**

P1410P HAND BED RAIL WITH POUCH



P1411 HAND BED RAIL WITH FLOOR SUPPORTS

- 20" wide handle has ¼" thick foam tubing to provide a secure grip
- Slides up to 43" under mattress with easy adjustments and has non-slip sleeves to keep it in place with adjustable strap to secure to frame
- Made of powder coated steel

SPECIFICATIONS

User Weight Limit: 300 lbs.
 Height Adj: 19¼" - 21¼"
 Leg Adj: 13¾" - 21"



Introducing the **NEW**

P1411P HAND BED RAIL WITH FLOOR SUPPORTS AND POUCH



P1409 BED RAIL POUCH

- Three pockets will fit remotes, cell phones, reading materials, etc.
- Fits almost any bed rail

To order, you may:

- Contact your Essential Medical Sales Rep.
- Call (800) 826-8423
- Visit www.essentialmedicalsupply.com