

Excel Kidz Pediatric

Features

- Adjustable telescoping handles for caregiver comfort
- Swing back desk length permanent arms
- Durable vinyl upholstery
- Padded and upholstered armrests and calf pads come standard
- Treaded, flat-free tires for safe riding
- Anti-tippers included
- Standard dual axle allows for easy adjustment of seat height to hemi level
- Flip-back arms

Accessories

Item Number	Description
MDS85181U	Oxygen Holder
MDS85183	I.V. Pole
MDS85190	Oxygen / I.V. Combo
MDS85180	5-in-1 I.V. / Oxygen / Anti-Theft
MDS85196K	Anti-Fold Device
MDS85197K	Overhead Anti-Theft Device
WCA806991	Tinnerman Legrest Locks

250 lb | 113 kg
weight capacity



MDS806140PEDE

Item Number	Seat W x D: Inches (cm)	HCPCS
MDS806140PEDE	14 x 12 (36 x 30)	E1238 + E0971 + K0195 or E0990

[Click here to view detailed product information and pricing.](#)

Call 800-397-5899 to order today!