

Face Shield

McKesson E1 Item Number: 1166883

MFR Number: 11000-SMC

EC Number: 16831200

Specifications:

Color: Clear Shield, Gray Headband

Dimensions: 1 x 1 x 9 inches

Fastening Type: Foam Headband with Adjustable Plastic Strap

Size: One Size Fits Most

Type: Anti-Fog

Usage: Disposable

Packaging: 100 per case



Features:

Plastic, adjustable belt-strap headband

Lightweight, anti-fog protection from splashing or spraying

Shield can be used in conjunction with a N95 mask, goggles, glasses or be worn alone without fogging up

Full face protection from fluid sprays and splashes

Optically clear with full peripheral vision

Solid foam provides comfort and protection from above

Shield made of optically clear polyester film

Head piece made of thick polyurethane

McKesson Medical-Surgical Inc.
Richmond, VA 23233

Questions? Call 1-800-777-4908

[mms.mckesson.com](https://www.mms.mckesson.com)

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information. The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.