






MIGHTY 400 POOL LIFT

PRODUCT WEIGHT: 150lbs [68kg]
MAX WEIGHT CAPACITY: 400lbs [181kg]
MIN ANCHOR SETBACK: 24"
MAX ANCHOR SETBACK: 44"
MAX WATER DRAFT: 12"
SEAT HEIGHT ADJUSTMENT: 4" IN 2" INCREMENTS
VERTICLE TRAVEL DISTANCE: 48"



***MEASUREMENT TAKEN WITH THE LIFT AT FULL LOAD**

***12" MAX DECK TO WATER**



UPRIGHT

EXTENDED

CONTINUOUS 360° ROTATION



9889 GARRYMORE LANE
MISSOULA, MT 59808
(406) 549-0769
FAX (406) 549-2602

NOTE:
*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

PRODUCT:		MIGHTY 400	
PART #:		F-MTY400	
SCALE:	1:28	DATE:	8/29/2019
DRAWN BY:	KR	REVISION:	02