

# CORFLO®

## NASOENTERIC

Nasogastric/Nasointestinal  
Feeding Tubes

with ENFit™ Connectors


CORFLO® Feeding Tubes from CORPAK MedSystems provide the largest selection of adult, pediatric and neonatal feeding tubes. The new CORFLO clear Y access port with ENFit™ connectors is designed to help prevent feeding tube misconnections and enhance safety.\*

### Your CLEAR connection to Enteral Safety!

- Clear material allows for visual inspection of connection with administration sets and for cleanliness before, during, and after use.
- Soft, pliable material and ergonomic, smooth shape – was designed with human factors engineering and customer input.
- Purple caps with easy-grip ridges and threaded caps provide secure connections. Purple helps identify the tube is for enteral use.

### FEATURES

- The ONLY feeding tube that can be used with the patented CORTRAK® Enteral Access System (EAS™)
- Exclusive CORPAK YELLOW tube design
- CORFLO® ULTRA's unique anti-clog exit port design is more than 300% larger than the internal diameter of the feeding tube, significantly reducing the occurrence of clogging
- CORFLO® CONTROLLER tubes are designed for maximum control during the insertion procedure, especially with endoscopic and fluoroscopic placements
- Polyurethane material provides larger inner lumen, which is excellent for long-term use
- Braided stylet design allows simple removal
- Rounded tungsten weights provide flexibility and increased comfort during placement
- Entire tube is radiopaque for easy visualization under X-ray
- Centimeter marks along the length of the tube to aid placement and check migration
- MRI SAFE (with stylet removed, MRI ≤ 3 TESLA)



Your  
**CLEAR**  
Connection  
to Enteral  
Safety

**CORPAK**  
MEDSYSTEMS

Ask for the **YELLOW** Tube...

\*ISO 80369-compliant.

U.S. PATENTS D750,757, D750,758, D750,759, AND OTHER PATENTS PENDING

#### CORFLO® ULTRA\*\*

Designed for insertions where anti-clogging and ease of insertion are critical.

#### CORFLO® CONTROLLER

Designed for maximum control during the insertion procedure, especially with endoscopic and fluoroscopic placements.

#### CORTRAK® 2

#### Enteral Access System (EAS™)

Created for real-time feeding tube placement confirmation.

CORFLO-ULTRA Lite NG	Non-wt. w/stylet
CORFLO-ULTRA Lite NG	Non-wt. w/o stylet
CORFLO-ULTRA Pedi NG	1 g pediatric wt. w/stylet
CORFLO-ULTRA NG	3 g wt. w/stylet
CORFLO-ULTRA NG	3 g wt. w/o stylet
CORFLO-ULTRA Pill NG	3 g wt. w/stylet
CORFLO-ULTRA 7 NG	7 g wt. w/stylet
CORFLO-CONTROLLER NG	3 g wt. w/stylet
CORFLO-CONTROLLER Pill NG	3 g wt. w/stylet
CORFLO-CONTROLLER 7 NG	7 g wt. w/stylet
CORTRAK 2 EAS NG	Non-wt. w/stylet

SIZE	5 FR	6 FR	8 FR	10 FR	12 FR
15" (38 cm)	42-1155*	42-1156*	42-1158*		
22" (56 cm)	42-9225*	40-9226	40-9228		
	42-1225*	40-1226	42-1226		
		42-1226*			
		40-8226			
36" (91 cm)	42-1365*	40-9366	40-9368TRAK2	40-9361TRAK2	40-1362
		40-1366	40-9368	40-9361	40-5362
		40-8366	40-1368	40-1361	
		40-7366	40-8368	40-7361	
			40-7368	40-5361	
43" (109 cm)			40-9438TRAK2	40-9431TRAK2	40-9432TRAK2
			40-9438	40-9431	40-9432
			40-1438	40-1431	40-1432
			40-7438	40-7431	40-7432
			40-5438	40-5431	40-5432
			40-6438	40-6431	40-4432
			40-4438	40-4431	40-0432
			40-3438	40-3431	
			40-2438	40-0431	
			40-0438	40-2431	
55" (140 cm)			40-9558TRAK2	40-9551TRAK2	
			40-9558	40-9551	
				40-3551	
				40-2551	

\* Sterile

\*\* All CORFLO® ULTRA feeding tubes have CORPAK's unique Anti-Clog exit port.

Boluses: approximate weights

All feeding tubes packaged 10/case



Innovation | Quality | Safety

#### TO ORDER Call Customer Service:

U.S.: Toll Free 1-800-323-6305 or 1-847-403-3400

FAX: 1-847-541-9526

For more product information, please visit  
[www.corpak.com](http://www.corpak.com)

CORPAK University Training  
& Education Portal:

[www.corpakU.com](http://www.corpakU.com)

