

1st Response™ Adult Manual Resuscitators



1st Response™ Features:

- Collapsible volume bag design for ease of storage.
- Low dead-space, swivel patient connector.
- Patient connector with integral swivel improves range of motion.
- 30 mm exhalation port.
- Tethered cap for clean, convenient storage of bag after use.
- Volume bag “full body” textures enhance security of grasp.
- Stroke volume of 850 ml.

PORTEX™

1st Response™ Adult Manual Resuscitators

ORDERING INFORMATION



PRODUCT CODE	DESCRIPTION	UNITS PER CASE
8500	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap Flow Diverter Reservoir Bag Adult Mask	12
8500P	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap Reservoir Bag PEEP Valve Adult Mask	12
8502	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap Reservoir Bag without Mask	12
8503	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap Flow Diverter 40" Flexible Tubing Reservoir Adult Mask	12
8503P	SureGrip™ Textured Resuscitator Bag with Tethered Dust Cap 40" (101 cm) Flexible Tubing Reservoir PEEP Valve Adult Mask	12
8505	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap 40" Flexible Tubing Reservoir without Mask	12
8506	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap Flow Diverter Expandable Tubing Reservoir Adult Mask	12
8506P	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap Expandable Tubing Reservoir PEEP Valve Adult Mask	12
8508	SureGrip™ Textured Adult Manual Resuscitator with Tethered Dust Cap Expandable Tubing Reservoir without Mask	12

All adult bags come with our standard 10ft (crush resistant) oxygen tubing.



Viewable Duckbill



Collapsible



Enhanced Texture



Exhalation Filter

Manual Resuscitator Accessories

ORDERING INFORMATION



PRODUCT CODE	DESCRIPTION	UNITS PER CASE
8501	Disposable PEEP Valve (Adjustable up to 20 cm H ₂ O) with 30 mm I.D. Connection	12
8518	Disposable Exhalation Filter with 30 mm I.D. Connection	40
008015	Resuscitation Valve Assembly with Filter	10



PEEP Valve

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.

160 Weymouth Street, Rockland, MA 02370
Phone: (800) 848-1757 Fax: (800) 621-2590

www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Please refer to the Instruction for Use supplied with the product for detailed instructions, warnings, and cautions. The Smiths Medical, Portex and 1st Response design marks are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © 2007 Smiths Medical family of companies. All rights reserved. RE435 10.2007



smiths medical