

1st Response™ PEDiplus® Manual Resuscitators



1st Response™ PEDiplus® Manual Resuscitator Features:

- The stroke volume is 600 ml.
- Low dead-space, swivel patient connector.
- Has manometer port so you have the ability to measure airway pressure.
- Longer tethered caps for clean, convenient storage of bag after use.
- Patient connector with integral swivel improves range of motion.
- Meets ASTM 920 and ISO 8382 standards.
- Optional Pressure Limiting Valve set at 40 cm H₂O (4kPa).

PORTEXTM

[Click here to view detailed product descriptions and prices.](#)

Call 800-397-5899 to order today!

1st Response™ PEDplus® Manual Resuscitator

ORDERING INFORMATION



PRODUCT CODE	DESCRIPTION	UNITS PER CASE
8515A	PEDplus® Large Pediatric Manual Resuscitator 40' Flex Tube Reservoir Small Adult Mask	9
8515C	PEDplus® Large Pediatric Manual Resuscitator 40' Flex Tube Reservoir Child Mask	9
8515P	PEDplus® Large Pediatric Manual Resuscitator 40' Flex Tube Reservoir PEEP Valve Small Adult and Child Mask	9
8515AB	PEDplus® Large Pediatric Manual Resuscitator Bag Reservoir Small Adult Mask	9
8515CB	PEDplus® Large Pediatric Manual Resuscitator Bag Reservoir Child Mask	9
8515PB	PEDplus® Large Pediatric Manual Resuscitator Bag Reservoir PEEP Valve Small Adult and Child Mask	9

All PEDplus bags come with our standard 10ft (crush resistant) oxygen tubing.



Viewable Duckbill



Collapsible



Exhalation Filter

Manual Resuscitator Accessories

ORDERING INFORMATION



PRODUCT CODE	DESCRIPTION	UNITS PER CASE
8501	Disposable PEEP Valve (Adjustable up to 20 cm H ₂ O) with 30 mm I.D. Connection	12
8518	Disposable Exhalation Filter with 30 mm I.D. Connection	40
008015	Resuscitation Valve Assembly with Filter	10
008201	Disposable Manometer	20



PEEP Valve



Disposable Manometer

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.

160 Weymouth Street, Rockland, MA 02370

Phone: (800) 848-1757 Fax: (800) 621-2590

www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

Please refer to the Instruction for Use supplied with the product for detailed instructions, warnings, and cautions. The Smiths Medical, Portex, 1st Response and PEDplus design marks are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © 2007 Smiths Medical family of companies. All rights reserved. RE436 - RE194201EN 2.2007



smiths medical