

CELOX™ GRANULES

Stops bleeding, saves lives

Celox™ Granules proven haemostatic control

- Significantly reduces blood loss¹
- Treats irregular shaped wounds and multiple injuries
- Easy to use



Contact us at www.celoxmedical.com




CELOX™ GRANULES

Stops bleeding, saves lives

Celox™ Granules

Celox™ Granules can stop even potentially life threatening bleeding fast. Celox™ has been proven to work on moderate to severe bleeding and has been safety testing to the intense class III CE mark standards. In an independent trial by the US Navy, when compared to Quikclot* and gauze, Celox™ was the only product to give 100% survival and was also the only product to give robust clotting with no rebleeding¹

How Celox™ works



- Chitosan in Celox™ absorbs blood fluid, swells and forms a gel.
- Celox™ electrostatically attracts red blood cells and forms a gel-like plug.
- Does not rely on the body's own clotting mechanism.
- No heat generated.

Easy to use

The Celox™ granules mould to the shape of the wound to get firm pressure on to the source of the bleeding even in complex shaped injuries.

- Works on hypothermic blood.
- Residual chitosan is naturally broken down and excreted by the body.

Use on:

- Arterial and venous bleeding.
- Life-threatening bleeding.
- Complex and deep wounds.

Product information

Celox™ granules are available in a 35g sachet for major injuries, even if the casualty has more than one bleeding wound. We also offer the 15g sachet for situations where every gram of weight saved is vital, still with enough material to treat a life-threatening injury.

Ordering Information

Celox™ Granules:
35g sachet: FG08830241
15g sachet: FG08830181

Medtrade Products Ltd.

Electra House, Crewe Business Park, Crewe, CW1 6GL UK
T: +44 (0)1270 500019 F: +44 (0)1270 500045

MT/13/009

Contact us at www.celoxmedical.com

1. Journal Academic emer med Jan 08, (Pgs 74-81). TM Trademark of Medtrade Products Ltd. *Trademark of Z-Medica Corporation.