


Smith & Nephew

PRIMAPORE[◇]

Conformable fabric island dressing with absorbent pad

PRIMAPORE is a soft, sterile dressing that can be used to cover a variety of wounds. It is ideal for use in treating wounds in curved areas of the body, including surgical wounds, general cuts, lacerations and abrasions. The highly absorbent pad provides an effective barrier, minimising fluid strike-through. PRIMAPORE also offers a low-adherent wound contact layer, reducing trauma to underlying tissue and ensuring pain is minimised on removal. The low-allergy acrylic adhesive provides secure fixation over the entire wound site.



Code	Product description (dressing; pad)	Items per unit
66003634	7.2cm x 5cm; pad: 4cm x 2.5cm	Box/50*
7135	8.3cm x 6cm; pad: 5.7cm x 3.5cm	Box/50
66000317	10cm x 8cm; pad: 5.3cm x 4.5cm	Box/20*
66000318	15cm x 8cm; pad: 11cm x 4cm	Box/20
66000319	20cm x 10cm; pad: 15.5cm x 5.3cm	Box/20
66000320	25cm x 10cm; pad: 20cm x 5.3cm	Box/20
66007140	35cm x 10cm; pad: 29.5cm x 5.3cm	Box/20*

[Click here to view detailed product information and pricing.](#)

Call 800-397-5899 to order today!