

PHILIPS

InnoSpire Essence

Nebulizer compressor system



Therapy you
can count on

innovation ✨ you



A **reliable** choice

Philips InnoSpire Essence nebulizer compressor system is designed for home-use and provides efficient and reliable aerosol therapy.

**There's always a way to
make life better.**

Reliable aerosol delivery you can trust



Compact design suitable for the home environment

InnoSpire Essence is for intermittent use (30 mins ON – 30 mins OFF) and is suitable for use in the home environment.



Compatible with commonly prescribed nebulizer medications

Capable of nebulizing a wide variety of medications that are routinely used with nebulizers to treat various respiratory conditions.



Treatment delivered in 6-8 minutes*

Combined with the SideStream Disposable Nebulizer, InnoSpire Essence is designed to provide fast medication delivery and short treatment times*.



Proven nebulizer technology

Essence is matched with our SideStream Disposable nebulizer to deliver a high-quality aerosol.

Find out more

www.philips.com/asthma



Proven technology for **fast and easy** treatment

SideStream has been designed to draw more air in, via an active venturi system, resulting in fast medication delivery and short treatment times. It also has a unique diamond jet design which helps to minimize wear from frequent use and aids consistent delivery of a high-quality aerosol. SideStream is suitable for use with most commonly prescribed respiratory medications.



SideStream has a unique venturi built into the nebulizer's top, which is designed to draw in extra air and shorten treatment times.



Innovative jet design producing consistently high-quality aerosol.



* Using 2.5ml salbutamol

[Click here to view detailed product information and pricing.](#)

Call 800-397-5899 to order today!



© 2019 Koninklijke Philips
N.V. All rights reserved.
Specifications are subject
to change without notice.
Respironics Inc. and
Respironics Respiratory Drug
Delivery (UK) Ltd are Philips
companies.

Respironics Respiratory
Drug Delivery (UK) Ltd
Chichester Business Park
City Fields Way
Tangmere
Chichester
PO20 2FT UK

PN 1139656
MCI 4108326
PN2805267
NP 12/08/19
InnoSpire Essence Consumer
Brochure, Gl Eng

www.philips.com/asthma