



## EMERGENCY OXYGEN UNITS & CPR ADMINISTRATION EQUIPMENT

[Emergency Oxygen Units](#)[Cylinders](#)[Regulators](#)[Masks](#)[Home](#)[AHA Emergency Oxygen Guidelines](#)[Contact LIFE®](#)[About Us](#)[Refill Information](#)[Maintenance](#)[Regulations](#)[Warranty](#)

### LIFE® SoftPac™

[Go To LIFE® StartSystem](#)[Go To LIFE® OxygenPac](#)

#### AED COMPANION O2 UNIT

Lightweight portable companion for your AED to provide Emergency Oxygen to a breathing victim before the onset of fibrillation, oxygen enriched CPR to a non-breathing victim, or continued supplemental oxygen after successful defibrillation.

- **Light Weight...Only 6 Pounds!**  
4 1/2" D x 7 1/2" W x 15" H
- Comfortable Handle & Shoulder Strap
- **Wall Mount Hook** for easy access
- **Constant reading supply gauge** through window (other units must be turned on to read gauge)
- **40 Minute Supply @ 6 Liters Per Minute**
- 250 Liter Aluminum "C" Cylinder
- **Cylinder is easily disposable/replaceable** or refillable at local gas distributors
- Simple External Two-step instructions  
**Inhalator:** administer supplemental oxygen for a breathing victim  
**Resuscitator:** oxygen enriched mouth-to-mask CPR for a non-breathing victim
- **5 Year Warranty**



*This unit could save your life.*



Make LIFE® SoftPac™ part of your first aid and safety program.

MODEL	FEATURES	Regulators
-------	----------	------------

LIFE-2-612	6&12 LPM-"NORM"&"HIGH"(2 settings) provides to recommended 100% inspired O2	
------------	---	--

LIFE-2-025	0 to 25 LPM Variable Flow (only for trained EMTs) 12 flow settings-0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25	
------------	--	--

LIFE-2-6FF	6 LPM (Liters Per Minute) Fixed-Flow	
------------	--------------------------------------	--

LIFE-2-101	<b>Replacement Cylinder Assembly</b> complete with cylinder, valve, gauge & O2 fill	
------------	--	--

*Uncompromisingly manufactured with quality materials, technology, and workmanship*



[LIFE® CPR Mask](#)  
invariably fits both ADULT and CHILD

#### "Knurled Knob" ON/OFF VALVE



Durable, Easy-to-Use, On-Off Valve  
(With Constant Reading Gauge)

[FDA Regulation](#) minimum capability is 6 LPM flow rate. LIFE® 612 models deliver both the minimum of 6 LPM and the AHA recommended 100% inspired oxygen at 12 LPM, in just two simple settings which read "NORM" & "HIGH". Only LIFE Corporation offers the simple 6&12 LPM flow rates.

#### American Heart Association Guidelines for CPR Oxygen

- ▶ "... give 100% inspired oxygen during Basic Life Support ... and Advanced Cardiac Life Support ... (BLS & ACLS) ... as soon as it becomes available."
- ▶ "... use supplemental oxygen...at a flow rate of 10-12 LPM ... to enable delivery of 100% inspired oxygen" (initially).

[www.LIFECorporation.com](http://www.LIFECorporation.com)

[Click here to view detailed product descriptions and prices.](#)

Call 800-397-5899 to order today!