

Compressor Nebulizer

Item#

NEB100

Comes with Disposable Adult Neb Kit and Filters

NEB100-R

Comes with Disposable and Reusable Adult Neb Kits and Filters

Neb kits include nebulizer bottle, t-piece, mouthpiece, 6in corrugated reservoir tube, and 7ft supply tube.

HCPCS# E0570 + A7003 + A7013

Item#

NEB100-BAG

Carry bag



Great value

**FAST
TREATMENT
TIMES**

Specifications

Weight 2.75 lbs

Input 120V/60Hz

Current 1.5A

Medication Capacity 6ml

Particle Size Respirable Fraction 0.5 to 5um 0.5um 80%

MMAD 2.0um

Average Nebulization Rate 0.3-0.4ml/min

Noise Level Below 60dBA

Compressor Pressure Range 30-45psi

Operation Pressure Range 12-14.5psi (0.85-1.1 bar)


Liter Flow Range 8-10LPM

Operation Temperature Range 50°F-104°F (10°C-40°C)

Operation Humidity Range 10 to 95% RH

Storage Temperature Range 14°F-149°F (-10°C-65°C)

Storage Humidity Range 10 to 95% RH



Comes in an attractive retail box

