

Personal Protective Gown



Knit polyester provides added comfort.

McKesson E1 Item Number: 1165820
MFR Number: GOWN01D50CP
AS400 Number: 58021100

Material: Polyester/Cotton (50/50)
Fabric Weight (g/m2): 190
Color: Yellow
Sterility: Non-sterile
Fluid resistant

Country of Origin: Honduras
Packaging: 12 Gowns/Case
Carton Dimensions (L x W x H):
 20" L x 12" W x 7.4" H

Size: One Size Fits Most
Caution: This product contains natural rubber latex which may cause allergic reactions.

MEASUREMENTS	in	cm
Height	44.5	113
Width	27.75	70.5
Length of Waist Tie	92	233.5
Length of Neck Tie	12.5	31.75
Length of Cuff Opening:	3	7.62
Size of Neck Opening:	9	22.86
Length of Sleeve:	22	55.5
Width of Sleeve:	11.5	29



McKesson Medical-Surgical Inc.
 Richmond, VA 23233

Questions? Call 1-800-777-4908

[mms.mckesson.com](https://www.mms.mckesson.com)

This product should not be used in a surgical setting or where significant exposure to bodily fluids or other hazardous fluids may be expected. This product should not be used in a clinical setting where Level 3 or 4 protection is warranted, or in any setting involving invasive procedures or where there is a high risk of contamination. This product should not be used in the presence of a high intensity heat source or flammable gas. This product has not been FDA cleared or approved. This product has been authorized by the FDA under the EUA for use as PPE in HC settings by HCP as they may help protect HCP and/or patients from the transfer of SARS-CoV-2 virus in low or minimal risk level situations to prevent the spread of COVID-19. This product is only authorized for the duration of the declaration.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information. The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.