


Hydraulic Lift

- Heavy gauge steel construction
- Six point spreader bar with 360° rotation
- Uses Lumex® six-point, four-point or two-point slings
- Foot pedal manually adjustable base
- Ergonomic placement of hydraulic pump handle facilitates ease of use
- Disassembles for easy storage
- Maximum weight capacity 400 lb
- Optional digital scale available
- Designed to meet the requirements of HCPCS code:E0630
- 1 year limited warranty

Specifications

Type.....	Hydraulic
Max Weight Capacity.....	400 lb
Floor Clearance.....	4.5"
Lifting Range Min Height.....	24" (75.5" max)
Max Height.....	77.7"
Overall Length.....	51"
Overall Height.....	55.1"
Overall Width.....	39.7"
Max Base Width.....	39.7"
Min Base Width.....	25.4"
Weight.....	90 lb
Warranty Frame.....	1 yr.

- Blue*
LF1030 1 ea
- Black/Silver Vein Hammertone*
LF1031 1 ea
LF1031FP with foot pedal adjustment 1 ea



Open



Closed



LF1030



LF1031FP