



A pool for everyone

Make Hippocampe chairs available at your pool!


- Access the changing rooms
- Cross the wet zone
- Access the showers, toilets, etc.
- Cross the foot bath
- Access the pool area
- Avoids annoying, tiring transfers from wheelchair to plastic chair when taking a shower.
- Access to the pool via an access ramp or steps (assistance provided).

A simple, innovative idea

- Water-resistant
- Easy to handle
- Seat made from hydrophobic foam (closed-cell polymer)
- Easy to maintain and disinfect - no damage posed to the wheelchair

«Swimming is the most popular sport among disabled persons, because it allows the majority of people with a moderate or severe handicap to use their motor functions under the condition of weightlessness»

Dealer :



Safety and Hygiene

The user goes into the shower on the Hippocampe chair before swimming. Thus the chair is regularly rinsed and disinfected, just like the other materials. It is completely rust-proof and is made from aluminium and plastic, which means that it is totally grease-free. The wheel bearings are made from self-lubricating erthacetal, which means there is no pollution. The chair seat is made from hydrophobic foam, which guarantees easy and effective rinsing.



A leisure area for everyone