

All-Terrain Knee Walker

HCPs:



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E

- Two 16" front wheels and two 12" rear wheels with enhanced traction effortlessly maneuver over gravel, dirt, grass, mulch, and uneven pavement
- Heavy-duty, extended-wear disc brake enables quick stopping for unparalleled safety
- Knee cushion with suspension is crafted for comfort to reduce potential pressure on shin and knee
- Foldable steering column allows for secure, compact storage and transport
- Complimentary water bottle and custom bottle holder make staying hydrated during travel easy

PRODUCTS

Description

| | |
|------|-----------------------------|
| 990X | All-Terrain Knee Walker, EA |
|------|-----------------------------|

SPECIFICATIONS

Description

| | |
|------------------------|------------------------|
| Front Wheel/Rear Wheel | Front: 16" Rear: 12" |
| Handle Height | 32.875" - 37.625" |
| Knee Pad Height | 18.375" - 23.25" |
| Warranty | Limited Lifetime |
| Weight | 35 lbs. |
| Weight Capacity | 350 lbs |

PARTS

Description

| | |
|------------|--|
| 9503H79002 | Hand Grips,790,1Pair, PR |
| 9503H79005 | Handle Bar Column,790,1ea, EA |
| 9503H79008 | Brake Cable Guide,790,1ea, EA |
| 9505B79006 | HandleBarHeightAdjLock, 790, 5B Ser#,1e, EA |
| 9503H79012 | Column To Fork Locking Bolt,790,1ea, EA |
| 9503H79007 | Steering Column,790,1ea, EA |
| 9503H79013 | Fork Bearings,790,1pair, PR |
| STDS1012 | Rear Wheel Bearing,11.5 x 32 x 9.5mm,1/e, EA |