

SecureSocks[™] For Fall Management

Model # - SOCKS-1B (black), SOCKS-1G (gray)
 SOCKS-1Y (yellow), SOCKS-1R (red) & SOCKS-1P (pink)

One Size Fits Most
 17" Length



Available in seven colors:
 Black, Gray, Yellow, Red & Pink

ALL AROUND NON-SKID PATTERN • SOFTER • LONGER

SecureSocks[™] were designed as a cost-effective complement to your fall management program. The socks provide an all around the sock non-skid pattern to assist in preventing falls with high fall risk patients/residents. Unlike competitive fall management socks, the SecureSock[™] was designed to maintain, as complete as possible, an all around non-skid pattern when placed on a patient/resident's foot. Compare the slippery gaps of competitive socks when placed on a foot with the SecureSock[™] all around non-skid pattern to see the safety difference.

SecureSocks[™] are made of ultra soft 80% cotton, 18% nylon & 2% Spandex for maximum comfort and warmth. The socks are longer than competitive socks to help prevent them from coming off a patient/resident's foot as a result of a shuffling walk.


Washing instructions:

Hand wash or machine wash on delicate cycle. Use warm water with a mild detergent. Do not use bleach. Tumble dry on low heat.

Warning!

Beware of wet floor surfaces. SecureSocks[™] were not designed to prevent falls when users are walking on wet or slippery floor surfaces.

Please see our complete line of Secure[®] fall & wandering prevention products at www.PadAlarm.com.




SecureSleeves[™]
 For Arms



SecureSleeves[™]
 For Legs



Secure[®] Hip Protectors - SHP
 (In five sizes & 3 colors)



Secure[®]
 Bed Pads-1 Set

Highest Quality...Lowest Prices...Guaranteed!

Personal Safety Corporation • Hiawatha, Iowa 52233 • 800.373.2873