



Standard



Re-Order #	Description	Size	Units
503100	Foot Elevator - Standard	4" Width	Each
503140	Foot Elevator - Standard	6" Width	Each

### Purpose

Designed for effective decubitus ulcer prevention and management. The Skil-Care Foot Elevator raises your patient's heel off the bed to eliminate pressure and also provides the necessary air circulation for accelerated healing of existing decubiti.

### Easy to Use

Wrap Foot Elevator around patient's ankle and adjust to a snug, comfortable fit. Secure hook-and-loop closure. There's no need to wrap a liner around the patient's ankle; the soft flannelette material will not irritate the patient's skin.

### Easy Care

This device may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

[Click here to view detailed product descriptions and prices.](#)

**Call 800-397-5899 to order today!**