


Visit us online
at Skil-Care.com

Skil-Care[™]
CORPORATION 29 Wells Avenue
Yonkers, NY 10701

914-963-2040
800-431-2972

SYNTHETIC SHEEPSKIN BED RAIL PADS

Re-Order #	Previous#	Size	Unit
□402010	9090	18" x 60"	Pair



Skil-Care's plush, high-pile synthetic sheepskin bed rail pads provide full protection by covering the patient side and top of the bed rail. Easy-to-use, heavy-duty snap secure the pad to the top of the rail. Quick-ties secure to the bottom. For total convenience, these pads freely raise and lower with the bed rails. Specially manufactured for the healthcare industry, our synthetic sheepskin is non-allergenic and non-irritating to sensitive skin.

LAUNDERING INSTRUCTIONS

These pads may be washed and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

Skil-Care[™]
CORPORATION

Advanced design for nursing and therapy.[™]

[Click here to view detailed product descriptions and prices.](#)

Call 800-397-5899 to order today!