





With the **Invacare**[®] **Walker**, providers, clinicians and consumers gain the stability they need at an economical price. Invacare Walkers offer a wide, deep frame with a large number of height adjustments and are offered in junior and adult sizes. The walker has a lower side brace for added stability and is stable, lightweight and easy to lift and maneuver. Walkers include PVC handgrips that provide comfort and long-lasting wear and anti-rattle silencers for quiet operation.

Invacare[®] Dual Blue-Release[™] Adult Walker Features

Lightweight, wide, deep aluminum frame Dual Blue-Release mechanisms provide both visual and audible "locked" cues Rigid dual side brace design Large number of height adjustments Anti-rattle "silencers" for quiet operation PVC handgrips for comfort and long lasting wear Compatible with most Invacare Walker Accessories (not to be used with Walker Tray or Walker Pouch)



Invacare® Single-Release Walker Features

Lightweight aluminum frame

One button to press and fold

Single button provides both visual and audible "locked" cues

Large range of height adjustments

PVC handgrips for comfort and long-lasting wear

Compatible with most Invacare Walker Accessories (not to be used with Walker Tray or Walker Pouch)



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

Invacare[®] I•Class[™] Paddle Walker Features

Composite lower side brace stiffens walker frame and adds stability Folding mechanisms are easy to reach from a seated position Lighter and easier to lift and maneuver PVS handgrips wear longer and provide greater comfort Walkers fold easily for transport or storage

HCPCS Code: E0135 Adult/Jr./Single E0143 Wheeled Adult/Jr.

The Invacare® Heavy-Duty Dual-Release Paddle Walker

is strong, stable and simple to use, helping individuals to walk with confidence and ease. This walker also is compatible with Invacare leg accessories.

Features

Compatible with model # 6375 bariatric 5" fixed wheels and leg extensions

The **Invacare® Bariatric Adult Walker** provides a wide and deep frame that can support an individual up to 700 lb. The dual-release buttons provide both visual and audible "locked" cues for the user. The durable aluminum frame provides maximum stability for consumers and the legs offer variable height adjustments.

Features

Compatible with Invacare model # 6372 bariatric 5" fixed wheels and leg extensions

HCPCS Code: E0148



2. Dual-Release Paddle button (model 6441-A)

3. Bariatric size

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.





The **Invacare® Two Step Walker** features a wide, deep frame with a large number of height adjustments and a lower side brace for added stability. It features two stages of grips to help the user ease from a sitting to standing position. The Invacare Two Step Walker is light, easy to lift and maneuver, and features foam hand grips for comfort and anti-rattle silencers for quiet operation.

Features

Combination single-release folding walker and sit to stand assistance device

Can also be used as a toilet safety frame

Adjustable height settings

Accessories currently not available

HCPCS Code: E0135



The **Invacare®** Adult Hemi Walker is designed for patients with limited or no dexterity in one arm or hand. This lightweight walker combines support and stability, yet is easy to lift and maneuver, even with one hand.

Features

Wider base provides more stability than a quad cane Anodized aluminum tubing, lightweight and maneuverable Compact size for easy storage and transport Extra handgrip provides ease in folding Slip-resistant rubber tips

HCPCS Code: E0135



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

Walker Specifications



Model #	Width Inside	Width Inside	Depth Opened	Depth	Height	Product	Weight	
	Hand Grips	Base Legs	At Base	Folded	Adjustment	Weight	Capacity	
6240-A	17"	21"	16"	4"	32"-39"	5 lb. 7 oz.	300 lb.	
6240-JR	17"	20"	16"	4"	25"-32"	5 lb. 5 oz.	300 lb.	
6240-5F	17"	21"	18.5"	4"	32"-39"	6 lb. 6 oz.	300 lb.	
6240-JR5F	17"	20"	18"	4"	25"-32"	6 lb. 12 oz.	300 lb.	
6281-A	17"	21"	17.5"	4.5"	30.375"-37.375"	6 lb. 2 oz.	300 lb.	
6281-JR	17"	21.5"	16.5"	4.5"	25.375"-32.375"	5 lb. 14 oz.	300 lb.	
6291-A	17"	21.5"	17.5"	4"	30.375"-37.375"	5 lb. 5 oz.	300 lb.	
6291-1	17"	21.5"	17.5"	4"	30.375"-37.375"	5 lb. 5 oz.	300 lb.	
6291-JR	17"	20.5"	16.5"	4"	25.375"-32.375"	5 lb. 1 oz.	300 lb.	
6291-3F	17"	21.5"	17.5"	4.5"	32"-38"	5 lb. 5oz.	300 lb.	
6291-JR3F	17"	20.25"	16.5"	4.5"	27"-33"	5 lb. 1oz.	300 lb.	
6291-5F**	17"	21.5"	17.5"	4.5"	33"-39"	5 lb. 6 oz.	300 lb.	
6291-JR5F**	17"	20.25"	16.5"	4.5"	28"-34"	5 lb. 6 oz.	300 lb.	
6291-HDA	20"	24"	17.5"	4"	30.375"-37.375"	6 lb. 11 oz.	500 lb.	
6441-A	20"	23.5"	16"	4.25"	32.25"-41.25"	7 lb.	700 lb.	
P421B	16.5"	21.5"	19"	5"	31"-36"	5 lb. 8 oz.	286 lb.	
6252	n/a	n/a	n/a	3"	29.5"-34.5"	2 lb. 14 oz.	250 lb.	



	Patient Height	Units	Limited
	Range*	Per Carton	Warranty
6240-A	5'4"-6'6"	4	5 Years
6240-JR	4'4"-5'7"	4	5 Years
6240-5F	5'4"-6'6"	4	5 Years
6240-JR5F	4'4"-5'7"	4	5 Years
6281-A	5'3"-6'4"	4	5 Years
6281-JR	4'4"-5'7""	4	5 Years
6291-A	5'3"-6'4"	4	5 Years
6291-1	5'3"-6'4"	1	5 Years
6291-JR	4'4"-5'7"	4	5 Years
6291-3F	5'4"-6'4"	4	5 Years
6291-JR3F	4'6"-5'6"	4	5 Years
6291-5F**	5'6"-6'6	4	5 Years
6291-JR5F**	4'4"-5'6"	4	5 Years
6291-HDA	5'3"-6'4"	2	5 Years
6441-A	5'3"-6'10"	2	5 Years
P421B	5'3"-6'4"	5	5 Years
6252	5'1"-5'10"	2	5 Years

* All patient height ranges are approximate. Each individual should always consult with his/her physician or therapist to determine proper

adjustment and usage.

** With 13" rear leg extensions. "F" designation in Model no. 3" or 5" fixed wheels with glide tips.



1. Walker Basket

Model # 6096 For model # 6240-A, 6240-JR, 6240-5F, 6240-JR5F walkers only

2. Walker Basket

Model # 6098 Compatible with all walkers except model # 6291-HDA, 6252, 6300-JRA, 6300-ATA

3. Walker Tray

Model # 6007 Compatible with all walkers except model # 6291-HDA, 6252, 6300-JRA, 6300-ATA, 6240, 6281 walkers

4. Walker Tray

Model # 6002 For model # 6240-A, 6240-JR, 6240-5F, 6240-JR5F walkers only

5. Walker Pouch

Model # 6015-4 Compatible with all walkers except model # 6291-HDA, 6252, 6300-JRA, 6300-ATA

6. Platform Attachment

Model # 6027 Compatible with all walkers except model # 6291-HDA, 6441-A, 6252 or 6300-JRA, 6300-ATA



Invacare[®] Court-Side Glide[™] Attachment Model # 6096

For model # 6240-A, 6240-JR, 6240-5F, 6240-JR5F walkers only



1. Glide Tips

Model # 6280 Packaged 24 per carton Durable construction for long wear





**Also available in retail packaging with a minimum purchase of four sets. Order as Model # 6270-4 and 6271-4.

- 1. 5" Swivel Wheels with 13" Cambered Invacare[®] SureGlide[™] Brakes Model # 6268
- 2. 5" Swivel Wheel Attachments with Rear Glide Tips Model # 6267
- 3. 3" Swivel Wheels with 13" Cambered Invacare[®] SureGlide[™] Brakes Model # 6272
- 4. **3" Swivel Wheel Attachments with Rear Glide Tips** Model # 6273
- 1. **3" Single Fixed-Wheel Attachments with Rear Glide Tips** Model # 6270**
- 5" Single Fixed Wheel Attachments with Rear Glide Tips Model # 6271** Use 13" rear leg extensions for maximum adjustability
- 3. Tall 3" Single Fixed-Wheel Attachments with Tall Rear Legs and Glide Tips

Model # 6278 with 13" rear leg extensions Adjusts to accommodate persons up to 6'6" tall

4. Bariatric Wheel Kit


Model # 6372 for use with 6441-A. 6375 for 6291-HDA not shown Dual 5" fixed wheels for added stability Large number of height adjustments



1. 10" Rear Wheel Brakes

Model # 6264 Use with 3" or 5" single fixed-wheel attachments

13" Rear Wheel Brakes Model # 6265 (not pictured) Use with 3" or 5" swivel wheel attachments or tall 3" single fixed-wheel attachments



- Tall Leg Extensions Model # 6275 Adjusts to accommodate persons up to 6' 6" tall
- 13" Rear Leg Extensions with Glide Tips Model # 6266 Use with swivel or fixed wheels

Invacare[®] Court-Side Glide[™] Specifications



Model # (Walker)	Patient Height Range	Units Per Carton
6291, 6240, 6281 w/ 5" fixed wheels	5' 6" - 6' 4"	6 Pair
6291, 6240, 6281 w/ 3" fixed wheels	5' 4" - 6' 4"	6 Pair
6291JR, 6240JR, 6281JR w/ 5" fixed wheels	4' 8" – 5' 6"	6 Pair
6291JR, 6240JR, 6281JR w/ 3" fixed wheels	4' 6" – 5' 6"	6 Pair
6300-JRA Adult to Junior WalkLite	4'11" – 5' 6"	6 Pair
6300-ATA Tall Adult WalkLite	5' 4" – 6' 0"	6 Pair

Accessory Specifications



Model # Patient Height Patient Height Limited Range Adult Range Junior Warranty min. – max. min. – max. Fixed-wheel 6270 31" – 37" 26" - 32" 5 Years 3" single fixed-wheels w/rear glide tips 33" – 39" 6271 28" - 33" 5 Years 5" single fixed-wheels w/rear glide tips 31" - 37" 6277 26" - 32" 5 Years tall 3" single fixed-wheels w/glide tips (not pictured) 6278 33" - 40" 31" - 35" 5 Years tall 3" single fixed-wheels w/tall rear leg ext. w/glide tips 6372* 33.5" - 43.5" N/A 5 Years double 5" fixed-wheels for 6441-A bariatric walker 6375* 33.5" - 43.5" N/A 5 Years double 5" fixed-wheels for 6291-HDA Swivel wheel 6267 37" - 40" 32" - 35" 5 Years 5" swivel wheels w/rear glide tips 37" – 40" 6268 32" - 35" 5 Years 5" swivel wheels w/13" cambered Invacare® SureGlide™ brakes 6269 34" - 40" 29" - 35' 5 Years 3" swivel wheels w/13" rear leg extensions w/glide tips 6272 34" - 40" 29" - 35" 5 Years 3" swivel wheels w/13" cambered Invacare SureGlide brakes 34" – 37" 6273 29" - 32" 5 Years 3" swivel wheels w/rear glide tips Leg extensions 6266 33" – 40" 29" - 35" 5 Years 13" rear leg extensions w/glide tips 33" – 40" 28" – 35" 5 Years 6275 tall front and rear leg extensions Rear brake attachments 6264 33" – 37" 28" - 32" 5 Years 10" rear wheel brakes 6265** 33" - 40" 29" - 35" 5 Years 13" rear wheel brakes

Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

©2012

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Lotus is a registered trademark of Inovo, Inc. Specifications are subject to change without notification. Form No. 06-015 Rev. 03/12 111253



**For use with 3" or 5" swivel wheels or tall 3" fixed wheels.