

JUXTA-LITE™

ANKLE-FOOT WRAP

The Juxta-Lite™ Ankle-Foot Wrap gives ankle and foot coverage. The wrap comes in two standard sizes and can fit almost any foot shape. The Breathe-O-Prene® material is cool, light-weight, and allows the skin to breathe. The SILVERtec™ fabric treatment can help eliminate odor and the growth of bacteria. Also, the Juxta-Lite™ Ankle-Foot wrap is machine washable.

SIZING AND PRICING

JUXTA-LITE™ ANKLE-FOOT WRAP	SMALL	MEDIUM	MSRP*
	38231017	38232017	\$36 [§]

* Prices may vary

ANKLE-FOOT WRAP ARCH CIRCUMFERENCE MEASUREMENTS		
	SMALL	MEDIUM
H	19-26cm	24-31cm



NOTE: The maximum foot arch circumference for the Juxta-Lite™ Ankle-Foot Wrap is 31 cm. If your patient exceeds this measurement, the Comfort EZ Single-Band Ankle-Foot Wrap™ is available for order (see below).

[§] \$10 pricing with the purchase of a Juxta-Lite™ Standard Legging



EACH PACKAGE INCLUDES

1 - Ankle-Foot Wrap

COMPONENTS

Nylon, Polyurethane, Spandex, Silver, Satin Acetate
This product does not contain Latex.

INDICATIONS

Venous Insufficiency, Venous Stasis Ulcers, Post Thrombotic Syndrome, Dependent Edema

CONTRAINDICATIONS

Severe Peripheral Arterial Disease, Decompensated Congestive Heart Failure, Septic Phlebitis, Phlegmasia Cerula Dolens, decreased or absent sensation in the leg, allergy to compression materials, Moderate Peripheral Arterial Disease, infection in the ankle and/or foot