

## Mayo-Instrument Stand

HCPCS:



- Casters provide additional mobility and convenience
- Removable stainless steel tray measures 19" x 13"
- Tray height adjusts from 29 ½" - 47" with lock

PRODUCTS	Description
13071	Instrument Stand,Mobile,5caster,1/cs

SPECIFICATIONS	Description
Base Dimensions	21" (W) x 21" (D)
Carton Shipping Weight	12 lbs.
Height	29.5" - 47"
Tray Dimensions	19" (W) x 13" (D)
Warranty	Limited Lifetime
Weight	10 lbs.

ACCESSORIES	Description
13035T	Tray for 13035, 13045,13071,1/ea