

Walkers

Decreased mobility doesn't have to mean the loss of independence. McKesson walkers are a convenient and affordable mobility solution for users who have difficulty with balance or standing, suffer arthritic pain, or experience fatigue or shortness of breath from walking. Walkers are commonly preferred for stability and balance. McKesson walkers feature lightweight construction and a push-button release mechanism for easy folding and storage. Choose McKesson walkers for durability, convenience, and comfort.

Features and Benefits

- Sturdy, 1-inch, aluminum A-frame construction is lightweight without sacrificing durability for patients up to 300 lbs
- Steel legs and side braces maximize durability and strength to support patients up to 500 lbs (applicable models only)
- Push-button, dual-release mechanisms for easy-folding using fingers, palms or side of hand to accommodate users' abilities and travel convenience
- Rubber-like tips made of polyurethane blend to provide better grip than standard walkers
- 5-inch front wheels and rear glide caps allow walker to slide easily and smoothly over most surfaces (only for walkers with wheels)
- Adjustable height and vinyl-contoured hand grips
- U-shaped frame provides greater clearance and the cross frame design offers improved stability
- Available in adult and junior sizes
- Not made with natural rubber latex



■ Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.









The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.


The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Walkers

ORDERING INFORMATION

STANDARD FOLDING WALKERS						
MFR #	DESCRIPTION	SIZE	INSIDE GRIP WIDTH	WEIGHT CAPACITY	HCPCS	QTY
146-10200-4 	Aluminum	Adult	17 in	350 lbs	E0135	4 ea/cs
146-10201-4 	Aluminum	Junior	17 in	350 lbs	E0135	4 ea/cs
146-10220-2 	Aluminum with Steel	Adult	19-1/2 in	500 lbs	E0148	2 ea/cs
FOLDING WALKERS WITH WHEELS						
146-10211-4 	Aluminum	Junior	17 in	350 lbs	E0143	4 ea/cs
146-10210-1 	Aluminum	Adult	17 in	350 lbs	E0143	1 ea/cs
146-10210-4 	Aluminum	Adult	17 in	350 lbs	E0143	4 ea/cs
146-10220-2WW 	Aluminum with Steel	Adult	20-1/2 Inch	500 lbs	E0149	2 ea/cs
FOLDING TRAVEL WALKERS WITH WHEELS						
146-RTL10263KDR 	Aluminum	Adult	17 in	300 lbs	E1043	1 ea/cs
ACCESSORIES						
146-RTL100014	Walker Ski Glides	N/A	N/A	N/A	N/A	1 pair, 4 pairs/cs

 **WARNING:** This product can expose you to bisphenol A (BPA), which is known to the State of California to cause birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

McKesson Medical-Surgical
 9954 Mayland Drive
 Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
 800.347.2456

Home Medical Equipment
 888.822.8111

Primary Care
 866.625.2679

McKesson Brands Support
 800.777.4908