# **M**SKESSON

## Canes

McKesson canes provide a convenient, affordable and lightweight mobility solution for those who experience pain, weakness or mild balance issues that make walking difficult. Quad canes are ideal for users who require more stability and support. Unlike most walkers, standard style canes and small base quad canes generally fit on an average depth stair. McKesson also offers convenient folding canes for ease of portability and storage. McKesson canes deliver ease of use, convenient height adjustability and reliable performance.

## **Features and Benefits**

- Standard style canes for enhanced maneuverability allowing use in narrow spaces such as staircases
- Aluminum construction ultra-lightweight solutions and steel options for users requiring added strength and durability
- Various ergonomic designs include round, T-shaped, and offset handle style options accommodate user's strength, hand size, and intended use
- Vinyl, foam/rubber, gel, silicone gel, and molded wooden grips for varying levels of comfort, durability, and user's grip strength
- Packaged for retail use
- Not made with natural rubber latex



### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

#### **Never Compromise on Care**<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2019 McKesson Medical-Surgical Inc.

# **M**SKESSON

### **ORDERING INFORMATION**

				STANDARD CANES				
MFR #	BASE TYPE	HANDLE TYPE	FRAME MATERIAL	HANDLE MATERIAL	HEIGHT	COLOR	HCPCS	QTY
146-10302-6*	Standard	Round	Aluminum	Vinyl Contoured	29.75 to 38.75 in	Chrome	E0100	1 ea, 6 ea/cs
146-RTL10306	Standard	Offset	Aluminum	Foam Rubber	30 to 39 in	Black	E0100	1 ea, 6 ea/cs
146-10303-6	Standard	Offset	Aluminum	Foam Rubber	30 to 39 in	Silver	E0100	1 ea, 6 ea/cs
146-RTL10307	Standard	Offset	Aluminum	Foam Rubber	30 to 39 in	Bronze	E0100	1 ea, 6 ea/cs
146-10305-6	Standard	Offset	Steel	Vinyl Contoured	29.75 to 37.75 in	Black	E0100	1 ea, 6 ea/cs
146-RTL10303PF	Standard	Offset	Aluminum	Foam Rubber	30 to 39 in	Pink Floral	E0100	1 ea, 6 ea/cs
146-RTL10372RC	Standard	Offset	Aluminum	Gel	30 to 39 in	Red Crackle	E0100	1 ea, 6 ea/cs
146-RTL10342	Standard	Round	Aluminum	Foam Rubber	28.75 to 37.75 in	Silver	E0100	1 ea, 6 ea/cs
STANDARD FOLDING CANES								
146-RTL10304	Standard	T-Shape	Aluminum	Wood Contoured	33 to 37 in	Black	E0100	1 ea, 12 ea/cs
146-RTL10304BZ	Standard	T-Shape	Aluminum	Wood Contoured	33 to 37 in	Bronze	E0100	1 ea, 12 ea/cs
146-RTL10304CE	Standard	T-Shape	Aluminum	Silicone Gel Grip	33 to 37 in	Glow in the Dark	E0100	1 ea, 6 ea/cs
				QUAD CANES				
146-10301F-4	Small	Offset	Steel	Foam Rubber	30 to 39 in	Chrome	E0105	1 ea, 4 ea/cs
146-RTL10310	Small	Offset	Steel	Foam Rubber	30 to 39 in	Black	E0105	1 ea, 4 ea/cs
146-10300-4 ^	Large	Offset	Steel	Vinyl Contoured	29 to 37.5 in	Chrome	E0105	1 ea, 4 ea/cs
				ACCESSORIES				
146-RTL10322BKB	Cane Tip	N/A	Rubber	N/A	3/4 in Diameter	Black	N/A	1 ea, 10 ea/cs
146-10374BK-8	Cane Tip	N/A	Rubber	N/A	7/8 in Diameter	Black	N/A	2 bx, 8 bx/cs
146-RTL10320BKB	Quad Cane Tip, Small Base	N/A	Rubber	N/A	1/2 in Diameter	Black	N/A	4 ea, 10 ea/cs
146-RTL10351	Four Point Cane Tip	N/A	Rubber	N/A	3/4 in Diameter	Black	N/A	1 ea, 10 ea/cs



\* **WARNING:** This product can expose you to diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information, go to www.P65Warnings.ca.gov. (**MFR# 10302-6**)

\*WARNING: This product can expose you to di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov. (MFR# 10300-4)

## **M**SKESSON

McKesson Medical-Surgical 9954 Mayland Drive Richmond, VA 23233 mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health 800.347.2456

Home Medical Equipment 888.822.8111

Primary Care 866.625.2679

McKesson Brands Support 800.777.4908

Rev. 08 03/19