

- Designed with reinforced tubing and/or wider bases to allow greater weight-bearing
- High weight-bearing capabilities

# GUARDIAN®



MDS86228XW – Double footplate and cane



### Bariatric Large Base Quad 500 lb (227 kg) Capacity

All steel construction with double-plated base allows for a 500 lb (227 kg) weight capacity. Contains a secure locking ring and dual-tubing insert for added strength and safety.



### Bariatric Small Base Quad 500 lb (227 kg) Capacity

All steel construction with double-plated base allows for a 500 lb (227 kg) weight capacity. Contains a secure locking ring and dual-tubing insert for added strength and safety.



### Single Point Bariatric Cane 500 lb (227 kg) Capacity

All steel construction with dual-tubing insert for added strength and security. Single point design gives additional balance support and weight-bearing capabilities.



| Item                               | Description  | Color | Weight Cap<br>Lbs (kg) | Adjustment<br>Inch (cm) | Footprint<br>Inch (cm) | Cane Weight<br>Lbs (kg) | Pkg  |
|------------------------------------|--------------|-------|------------------------|-------------------------|------------------------|-------------------------|------|
| <b>Bariatric Quad Canes</b>        |              |       |                        |                         |                        |                         |      |
| MDS86228XW                         | Large Base   | Black | 500 (227)              | 29-38 (74-96.5)         | 8 x 12 (20 x 30)       | 4.6 (2)                 | Each |
| MDS86222XW                         | Small Base   | Black | 500 (227)              | 29-38 (74-96.5)         | 6 x 8 (15 x 20)        | 3.6 (2)                 | Each |
| <b>Single Point Bariatric Cane</b> |              |       |                        |                         |                        |                         |      |
| MDS86420XW                         | Single Point | Black | 500 (227)              | 29-38 (74-96.5)         | N/A                    | 3.0 (1)                 | Each |

## Cane Tips and Grip

- Tips are reinforced with metal inserts for longer wear

| Item Description                        | Size-Inches (cm) | Color | Pkg      |
|---|------------------|-------|----------|
| <b>Tips</b>                             |                  |       |          |
| G30561-PC Tip                           | 5/8 (1.6)        | Black | 8 pr/cs  |
| G30572-PC Tip                           | 3/4 (2.0)        | Gray  | 8 pr/cs  |
| G30581-PC Tip                           | 7/8 (2.0)        | Black | 8 pr/cs  |
| G02162 Tip, Quad Cane                   | 1/2 (1.3)        | Gray  | 8 set/cs |
| G02172 Tip, Quad Cane                   | 5/8 (1.6)        | Gray  | 8 set/cs |
| MDS80212T Tip, Replacement for MDS80212 | 7/8 (2.0)        | Black | 2 ea/pr  |
| MDS80217T Tip, Replacement for MDS80217 | 1 (2.5)          | Black | 2 ea/pr  |
| MDS86424W Tip, Quad Cane                | 1/2 (1.3)        | Black | 6 pr/cs  |
| MDS86425W Tip, Quad Cane                | 5/8 (1.6)        | Black | 6 pr/cs  |
| MDS86426W Tip, Quad Cane                | 3/4 (2.0)        | Black | 6 pr/cs  |
| <b>Grip</b>                             |                  |       |          |
| MDS86427W Grip, for 420,222,228         | 4 (10.2)         | Black | 6 pr/cs  |

