

# MONOJECT™ BLUNT FILL NEEDLE



The Monoject™ blunt fill needle is a reliable solution for medication preparation. Blunt fill needles are suitable for accessing I.V. bags, drawing-up and reconstituting medicine. The Monoject™ blunt fill needle is available with and without a 5 micron filter.

## FEATURES & BENEFITS

### Standard Luer Hub

Fits all standard luer lock syringes

### 45° Bevel Needle

Helps reduce needle-stick injuries

### 5 Micron Filter

Effectively filters contaminant particles for safe ampule medication preparation

### Red Identification Cap

For quick and easy identification



Attach the blunt fill cannula to a luer lock syringe to draw-up medication

# MONOJECT™

## BLUNT FILL NEEDLE

### PRODUCT ORDERING INFORMATION

Code	Needle Gauge and Length	Box/Case Qty
11811022	Bluntfill 18 G x 1-1/2" (1.27 mm x 3.8 cm)	100/1000
11811022F	Bluntfil with filter 18 G x 1-1/2" (1.27 mm x 3.8 cm)	100/1000

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2015 Medtronic. All rights reserved. Medtronic, Medtronic logo and *Further, Together* are trademarks of Medtronic.

™\* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

– EU-15-428458 – 10/2015

**Medtronic**