

Multi-Purpose Containers

Rotor Opening Lids



8970

Hinged Opening Lids



8990SA

Rotor & Hinged Opening Lids



8527R

Multi-Purpose Sharps Containers

- Adjustable rotor opening or hinged lids accommodate a variety of sharps sizes and provide temporary and permanent closure.
- Easy, simple assembly.
- Containers lock for final disposal.
- Nestable containers save valuable storage space.
- Available holders and brackets secure containers wherever needed.
- Plastics and pigments are safe to autoclave or incinerate.
- For use on countertops or mounted on the wall.
- Economical horizontal drop maximizes container volume.
- Large horizontal openings accommodate larger syringes.

Code	Description	Dimensions	Ship Case
Multi-Purpose Containers with Rotor Opening Lid			
8920SA	2 Quart, Red	4.58"H x 6.285"D x 4.75"W	60
8950SA	5 Quart, Round, Red	6.75"H x 8.75"D	40
8970	2 Gallon, Red	10"H x 7.25"D x 10.5"W	20
8979C	2 Gallon, Transparent Red	10"H x 7.25"D x 10.5"W	20
Multi-Purpose Containers with Hinged Lid			
8990SA	2 Gallon, Red	10"H x 7.25"D x 10.5"W	20
Multi-Purpose Containers with Rotor & Hinged Opening Lid			
8527R	3 Gallon, Transparent Red	13.75"H x 6"D x 13.75"W	10