





## Upper Arm Blood Pressure Monitor Comparison Chart

							
Model Number	3 Series BP7100	5 Series BP742N	7 Series BP760N	7 Series Wireless BP761N	10 Series BP7450	Evolv BP7000	Complete BP7900
Advanced Accuracy	✓	✓	✓	✓	✓	✓	✓
Hypertension Indicator					✓		
Bluetooth®				✓	✓	✓	✓
TruRead					✓		
Extra-Large Digits w/Backlight					✓		
One-Touch Use	✓	✓				✓	
Guest Mode		✓	✓	✓	✓		
Two User Mode		✓	✓	Unlimited users in Connect app	Unlimited users in Connect app	Unlimited users in Connect app	Unlimited users in Connect app
Body Movement Detection	✓	✓	✓	✓	✓	✓	
Advanced Averaging		✓	✓	✓	✓		
Irregular Heartbeat Detector	✓	✓	✓	✓	✓	✓	✓
Memory per User	14	50	60	60	100	Unlimited users in Connect app	Unlimited users in Connect app
Morning Hypertension Indicator							
Morning Average/Week							
EKG							✓
Side-by-Side Comparison Display				✓	✓		
Cuff	Wide-Range D-Ring 9"-17"	Wide-Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Easy-Wrap ComFit 9"-17"	Easy-Wrap ComFit 9"-17"	Easy-Wrap ComFit 9"-17"	Wide-Range D-Ring 9"-17"
Replacement Cuff	CD-WR17	CD-WR17	CFX-WR17	CFX-WR17	HEM-FL31-B	N/A	HEM-RML31-B
Optional Cuff	CD-CS9 (7"-9")	CD-CS9 (7"-9")					HEM-CS24-B (7"-9")
Storage Case							
Power Options	4 AA Not Included	4 AA Included / AC Adapter Not Included	4 AA Included / AC Adapter Not Included	4 AA Included / AC Adapter Not Included	4 AA Included / AC Adapter Included	4 AAA Included / Not compatible with AC Adaptor	4 AA Included / Not compatible with AC Adaptor