

CENTURION®



# SORBAVIEW® SHIELD



*"You can't pull out the catheter no matter how hard you tug.  
The SHIELD is a huge cost savings and works great."*

- S.S., PICC RN\*

# SORBAVIEW® SHIELD


**A one-step catheter securement system that combines the features and benefits of a SorbaView Dressing with SHIELD Technology.**

Centurion is inspired by real-world insight in the development of products that improve patient care. We work directly with doctors and nurses to understand their needs and produce products that meet their specific performance requirements. SorbaView SHIELD is a direct result of these relationships.

The SHIELD is a combination of our industry-leading SorbaView Dressing line and SHIELD Technology. **Our dressings go on easier and wear longer** than other dressings, resulting in fewer and faster changes. SorbaView's patented design wicks away excess moisture and "breathes" like skin, while protecting patients by keeping out bacteria. This helps lower the risk of catheter-related bloodstream infections (CRBSI) associated with central lines.

9

## VECTOR FORCE TUG



## SHIELD TECHNOLOGY

**Manufactured stabilization device built into a SorbaView Dressing that can withstand a 9-Vector Force Tug, preventing catheter movement and dislodgement.**

## SHIELD Technology guards the catheter from aggressive movements

and eliminates the need for a secondary securement device. With the addition of SHIELD Technology, dressings stay on and catheters stay put. By keeping the catheter in place, the risk of complications such as dislodgement, phlebitis, thrombosis and infiltration is dramatically reduced in IVs.

With the SHIELD, Centurion has created **the ultimate catheter securement solution.** Our premium-quality window dressing combined with SHIELD Technology results in a manufactured stabilization device that prevents catheter movement and dislodgement, while creating an incredibly strong barrier that protects the insertion site. Marrying these defensive features together into a one-step catheter securement system results in better patient care and lower costs.



*"We switched to the SorbaView in January. It is performing extremely well. The securement feature is great. We are also finding that the dressing adherances last the full 7 days, decreasing our premature dressing changes. Our staff loves the dressings ... there are NO complaints!"*  
- Donna S., IV Therapy Coordinator\*

- SorbaView SHIELD is founded on the solid 17-year successful history of the SorbaView product line.
- SorbaView Evidence Based Nursing data has been featured in several studies and published articles in Nursing Journals.<sup>†</sup>
- The SHIELD is manufactured in accordance with FDA regulations and meets the Infusion Nurses Society (INS) recommendations for manufactured stabilization devices.<sup>††</sup>

*"I tried SorbaView SHIELD, and within one week found my skin less irritated and itchy. The SorbaView SHIELD dressing made [a] huge impact on my skin! It was a night-to-day difference."*

*- Molly, leukemia patient\**



## SUPERIOR PERFORMANCE BY DESIGN

### DISCOVER THE BENEFITS OF THE CENTURION SORBAVIEW® SHIELD PRODUCT LINE

#### CAREGIVER BENEFITS

- SHIELD Technology provides effective catheter securement as part of the normal dressing procedure – no extra steps or extra time needed
- Easy to apply – won't stick to itself during application
- Combining two products into one makes application faster and easier than with two separate products
- IV lines stay put, even when patients are transferred between hospital departments, saving caregivers time by reducing the number of unscheduled restarts
- 7-day wear time means fewer dressing changes on central lines
- Available in multiple sizes to fit any vascular access site

#### HOSPITAL BENEFITS

- Establishes a standard of care and eliminates inconsistent methods of securement
- Wear time of 7 days meets the demands of hospital protocols for central lines and results in fewer dressing changes and reduced labor costs
- Provides proven bacterial protection, lowering the risk of contracting a hospital-acquired infection (HAI)
- Prevents catheter movement, so the risk of an IV complication, such as dislodgement, phlebitis, thrombosis or infiltration – is dramatically decreased
- Using one product instead of two lowers costs
- Available in individual sterile packages or in custom procedure trays for ordering convenience

#### PATIENT BENEFITS

- Prevents catheter movement and dislodgement, increasing patient comfort
- Latex free, non-sensitizing adhesive makes the SHIELD comfortable for patients to wear
- Provides protection against bacteria, lowering the risk of contracting an HAI
- Offers a positive patient experience by lowering the number of times the site is disturbed

*"I recently had an opportunity to try the SorbaView SHIELD on a patient that was having a problem with tape/Tegaderm™ skin irritation. We applied the product, and the patient, who was a nurse, loved it and her skin tolerated it well. She has been using it for close to a month now with no issues. The patient also commented that it was much more comfortable and softer feeling than the Tegaderm."*

*- T.B., Infusion Support Nurse\**

# UNIVERSAL APPLICATIONS



Arterial (SV254UDT)



Subclavian (SV254AFXT)



Jugular (SV353AFXT)



PICC (SV37UDT)



Pediatric Subclavian (SV254UDT)



Femoral (SV37UDT)

*"I can't tell you how happy I am with these dressings! Thanks."*

*- Chris, RN\**



Pediatric PICC (SV254UDT)



Implanted Port (SV254UDT)




Epidural (SV37UDT)

## NEW SIZE DESIGNED FOR PERIPHERAL APPLICATIONS

WORKS WITH MULTIPLE PERIPHERAL CATHETER TYPES



Peripheral (SV233UDT)



**LARGE**  
SV37UDT

**MEDIUM**  
SV353UDT

SV353AFXT

*With Adhesive-free Zone in Window*

**SMALL**  
SV254UDT

SV254AFXT

*With Adhesive-free Zone in Window*

**PERIPHERAL**  
SV233UDT

\*Quotes that appear throughout this brochure were unsolicited end-user and patient testimonials. Testimonial statements are anecdotal observations not necessarily supported by clinical research performed by Centurion Medical Products.

†References on file.

††Reference: Infusion Nursing Standards of Practice , Vol. 29, No. 1S

MADE IN THE U.S.A.

SorbaView® SHIELD and SHIELD Technology are the proprietary property of Centurion Medical Products.

U.S. Patent Nos. 6,841,715 and 7,294,752

**CENTURION**  
MEDICAL PRODUCTS

100 CENTURION WAY | P.O. BOX 510 | WILLIAMSTON | MICHIGAN | 48895  
517.546.5400 | 800.248.4058 | fax 517.546.9388 | www.centurionmp.com | www.sorbaviewshield.com