


SHILEY™ ACCESSORIES



- 8A. Inner Cannula, Disposable (DIC) with integral 15mm snap-lock connector, individually sterile packaged for patient protection, required for use with DCT, DFEN, PERC, DCFS and DCFN.
- 8B. Spare Inner Cannula (SIC) with red, color-coded twist-lock connector, for temporary use during cleaning or removal of the inner cannula of Shiley™ Tracheostomy Tubes, may be used with LPC, FEN, CFS and CFN.
- 8C. Inner Cannula, Fenestrated, Disposable (DICFEN) is for use with DCFN and DFEN tracheostomy tubes.
- 8D. White 15mm cap (CAP) attaches to the Fenestrated Inner Cannula with green 15mm connector to permit breathing through the fenestration and upper airway.
- 8E. Decannulation Plug (DCP) closes the proximal end of the FEN and CFN tubes to permit breathing through the fenestration and upper airway.
- 8F. Decannulation Plug, Disposable (DDCP) closes the proximal end of the DCFS, DCFN and DFEN to permit breathing through the fenestration and upper airway.

SIZE	PRODUCT DESIGNATION	FOR USE WITH				QUANTITY/ BOX	
4	4DIC	4DCT	4DFEN	4DCFS	4DCFN	10	DIC
6	6DIC	6DCT	6DFEN	6DCFS	6DCFN	10	
8	8DIC	8DCT	8DFEN	8DCFS	8DCFN	10	
10	10DIC	10DCT	10DFEN	10DCFS	10DCFN	10	
4	4SIC	4LPC	4FEN	4CFN	4CFS	10	SIC
6	6SIC	6LPC	6FEN	6CFN	6CFS	10	
8	8SIC	8LPC	8FEN	8CFN	8CFS	10	
10	10SIC	10LPC	10FEN	10CFN	10CFS	10	
4	4DICFEN	4DFEN	4DCFN			10	DICFEN
6	6DICFEN	6DFEN	6DCFN			10	
8	8DICFEN	8DFEN	8DCFN			10	
10	10DICFEN	10DFEN	10DCFN			10	
4	4DCP	4FEN	4CFN			10	DCP
6	6DCP	6FEN	6CFN			10	
8	8DCP	8FEN	8CFN			10	
10	10DCP	10FEN	10CFN			10	
	DDCP	Any size DCFS, DCFN and DFEN tracheostomy tube.				10	DDCP
	CAP	Any size FEN and CFN tracheostomy tube.				10	CAP