


# NITRO HD

## 10266HD-BK



- 10. Wheel Lock Assembly
  - Plastic Left- 1026641
  - Right- 1026642
- 11. Rear Wheel Assy. Right- 1026612-R  
Left- 1026612-L
- 12. Seat Upholstery & Hdwr.- 10266HD08
- 13. Strap- 1026609-FSG
- 14. Folding Locking Latch- 1026644
- 15. Back Height Adj. Lever- 1026607-R
- 16. Rear Axle Bearing (pr)- 1026613
- 17. Rear Hub Cap (pr)- 1026619
- 18. Rear Axle- 1026631
- 19. Front Axle Bearing (pr)- 1026618-AN
- 20. Front Hub Cap- 1026614-AN
- 21. Dust Cap (pr)- 1026620
- 22. Fork Assembly
  - Left- 1026622-AN
  - Right- 1026621-AN
- 23. Fork Stem Bearings (pr)- 1026623
- 24. Cross Frame Assist (pr)- 10266HD32
- 25. Front End Cap (pr)- 1026628
- 26. Rear End Cap (pr)- 1026629
- 27. Shock Absorber- 10266HD31
- 28. Step Tube Left- 1026636  
Right- 1026637
- 29. Brake Cable Guide- 1026639
- 30. Cane Holder- 892
- 31. Cross Brace Screw & Insert- 1026643
- 32. Front Axle- 1026645
- 33. Cross Frame- 10266HD-33  
Cross Frame Bolt- 1026646

- 1. Handgrip Right- 1026603  
Left- 1026602
- 2. Hand Brake Assy. Left- 10266HD01-L  
Right- 10266HD01-R  
Cable only- 10266HD02 (51" in length)
- 3. Handle Bar- 1026630N
- 4. Backrest Assy.- 10266HD04
- 5. Backrest Knob & Bolt (pr)- 1026605-G
- 6. Safety Lock- 1026638
- 7. Bag Mount & Screw (pr)-1026611
- 8. Bag- 10266HD10
- 9. Front Wheel & Fork Assembly
  - Left- 1026625-AN
  - Right- 1026624-AN
  - Front Wheel- 1026617-AN