




First Class School Chair

HCPCS: E1399



FIGURE A



FC 2000

- Seat – height and depth adjustable solid seat with a 15° anterior or posterior tilt
- Armrests – provide lateral limitations as well as a stable resting surface for arms
- Flip-up armrests (Figure A) allow easy transfers and they are height adjustable
- The tray accessory mounts onto the armrest
- Padded Hip Belt – 2-point hip belt provides pelvic support and maintains proper hip positioning

PRODUCTS	Description
FC 2000N	FirstClass SchoolChair Sm NewStyle, 1/bx
FC 4000N	22rstClass SchoolChair Lg NewStyle, 1/bx

SPECIFICATIONS	Description
Seat Width (FC 2000N FC 4000N)	12" 16"
Seat Depth (FC 2000N FC 4000N)	9.5"-14.5" 13"-17"
Seat to Floor Height (FC 2000N FC4000N)	10-14" 15 - 23"
Seat to Footplate (FC 2000N FC4000N)	9.5" - 12.5" 9.5" 17.5"
Backrest Height (FC 2000N FC 4000N)	15" 17"
Armrest Length (FC 2000N FC 4000N)	10" 10"
Seat-to-Armrest Height (FC 2000N FC 4000N)	5" - 10" 4" - 13.5"
Distance Between Armrests (FC 2000N FC 4000N)	13.5" 17.5"
Hip Guides (FC 2000N FC 4000N)	3" (H) x 10" (D) 3" (H) x 10" (D)
Distance Between Hip Guides (FC 2000N FC 4000N)	8" - 15" 10.5" - 19"
Lateral (FC 2000N FC 4000N)	4" (H) x 5" (D) 4" (H) x 5" (D)
Seat to bottom of Lateral Pads (FC 2000N FC 4000N)	3.5" - 8.25" 3" - 11.5"
Distance Between Laterals (FC 2000N FC 4000N)	6" - 13" 10" - 17"
Headrest Added Height (FC 2000N FC 4000N)	4.5" - 8" 4.5" - 8"
Product Weight (FC 2000N FC 4000N)	13 lbs. 18 lbs.
Shipping Carton Weight	16.5 lbs. 22.5 lbs.
Weight Capacity (FC 2000N FC 4000N)	100 lbs. 200 lbs.

ACCESSORIES	Description
FC 2028	Footrests/ FC 2000N, 1pr/bx
FC 4028	Footrests/ FC 4000N, 1pr/bx
FC 2024	Tray/ FC 2000N
FC 4024	Tray/ FC 4000N, 1/bx
FC 8000N	Headrest/ FC 2000N, FC 4000N, 1/bx
FC 2029N	Legs with Casters/ FC 2000N, 1set/bx
FC 4029N	Legs with Casters/ FC 4000N, 1set/bx
FC 2030N	Anti-Tippers/ FC 2000N, 1/bx

ACCESSORIES	Description
FC 4030N	Anti-Tippers/ FC 4000N, 1/bx
FC 8027N	Hip Guides/ FC 2000N, FC 4000N, 1pr/bx
FC 8025N	Support Kit/ Small, 1set/bx
FC 8025 LN	Support Kit/ Large, 1set/bx

PARTS	Description
FC 2002N	
FC 4002N	
FC 8018N	1/pr
FC 8060	3/ea
FC ARKNOB	1/ea
FC 8019N	1/ea
FC 2001N	1/ea
FC 4001N	1/ea
FC 2011N	
FC 4011N	
FC 2003N	
FC 4003N	
FC0002N	1/ea
FC 0003N	1/ea
FC 0014N	1/ea
FC KNOBN	4/ea