


## PVC POOL ACCESS CHAIR - 21 INCH WIDE MESH SEAT





## **SPECIFICATIONS:**

THE PVC POOL ACCESS CHAIR WITH 21 INCH WIDE MESH SEAT SHALL BE CONSTRUCTED OF HEALTH GRADE PVC PIPE AND FITTINGS THAT WON'T RUST OR STAIN. IT WILL BE RATED FOR 350 lbs [159 kg]. THE CHAIR SHALL COME EQUIPPED WITH AN EXTENDEDABLE FOOTREST. THE MESH SEAT SHALL BE A MINIMUM OF 21" WIDE AND UPHOLSTERED FROM MOLD AND MILDEW RESISTANT MATERIAL THAT IS SUITABLE FOR IMMERSION IN A POOL OR SPA. THE CHAIR SHALL COME EQUIPPED WITH DUAL LOCKING BRAKES AND SWING AWAY ARMRESTS TO ALLOW FOR EASY TRANSFERS. THE REAR WHEELS SHALL BE A NON-SKID STYLE WITH USER HAND RINGS. THE CASTERS SHALL BE \$\Infty\$ 3" MRI RATED SWIVAL TYPE.





9889 GARRYMORE LANE MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602

24" NON-SKID WHEELS

NOTE: 
\*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

PRODUCT:	F-520WSPM		
DWG#:	1 OF 1		
SCALE:	1:20	DATE: 4/12/2017	
DRAWN BY:	EF	REVISION: A	