


PVC POOL ACCESS CHAIR -24 INCH WIDE MESH SEAT





SPECIFICATIONS:

THE PVC POOL ACCESS CHAIR WITH 24 INCH WIDE MESH SEAT SHALL BE CONSTRUCTED OF HEALTH CARE GRADE PVC PIPE AND FITTING THAT WON'T RUST OR STAIN. IT WILL BE RATED FOR 425 lbs [193 kg]. THE CHAIR SHALL COME EQUIPPED WITH AN EXTENDEDABLE FOOTREST. THE MESH SEAT SHALL BE A MINIMUM OF 24" WIDE AND UPHOLSTERED FROOM MOLD AND MILDEW RESISTANT MATERIAL THAT IS SUITABLE FOR IMMERSION IN A POOL OR SPA. THE CHAIR SHALL COME EQUIPPED WITH DUAL LOCKING BRAKES AND SWING AWAY ARMRESTS TO ALLOW FOR EASY TRANSFERS. THE REAR WHEELS SHALL BE A NON-SKID STYLE WITH USER HAND RINGS. THE CASTERS SHALL BE \emptyset 3" MRI RATED SWIVAL TYPE.



24" NON-SKID WHEELS

Ø 3" MRI SWIVAL CASTERS



9889 GARRYMORE LANE MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602

*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

PRODUCT: F-520XWSPM		
DWG#: 1 OF 1		
SCALE:	1:20	DATE: 4/12/2017
DRAWN BY:	EF	REVISION: