

Cuffed Flex D.I.C® Tracheostomy Tubes

PORTEX®

The Portex® flexible disposable inner cannula tracheostomy tubes are constructed of siliconized PVC, providing a softer and more pliable tube which is intended to reduce risk of tracheal trauma. The built-in 15mm connector allows ventilator circuit attachment with or without inner cannula in place.



ORDERING INFORMATION

Box quantity: 1

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

PRODUCT BENEFITS

- Cuffed Flex D.I.C.® tubes come with a Soft-Seal® cuff, providing a high volume, low-pressure seal that reduces the potential for patient trauma.
- Tubes are available in five color-coded adult sizes. Color-coding throughout the D.I.C.® line enables caregivers to quickly identify and size tubes and related inner cannula.
- Obturator features special clip design to minimise tube tip movement during insertion.

PRODUCT SPECIFICATIONS

Cuffed Flex D.I.C.® Tracheostomy Tubes

Reorder Code	Size	ID (mm)	OD (mm)	Length (mm)	Colour Code
--------------	------	---------	---------	-------------	-------------

505060	6.0 mm	6.0 mm	8.2 mm	64.0 mm	Orange
--------	--------	--------	--------	---------	--------

505070	7.0 mm	7.0 mm	9.6 mm	70.0 mm	Green
--------	--------	--------	--------	---------	-------

505080	8.0 mm	8.0 mm	10.9 mm	73.0 mm	White
--------	--------	--------	---------	---------	-------

505090	9.0 mm	9.0 mm	12.3 mm	79.0 mm	Blue
--------	--------	--------	---------	---------	------

505100	10.0 mm	10.0 mm	13.7 mm	79.0 mm	Yellow
--------	---------	---------	---------	---------	--------

Notes: Flex D.I.C.® tracheostomy tubes come with an obturator, Velcro® brand tube holder, and one colour-coded inner cannula

Replacement Inner Cannula

Reorder Code Tube ID (mm) units / case

526060 5.0mm 60

526070 6.0mm 60

526080 7.0mm 60

526090 8.0mm 60

526100 9.0mm 60