

## Uncuffed Flex D.I.C® Tracheostomy Tubes

PORTEX®

The Portex® flexible disposable inner cannula tracheostomy tubes are constructed of siliconized polyvinyl chloride, intended to provide a softer and more pliable tube intended to reduce risk of tracheal trauma. Built-in 15 mm ISO termination for ventilator circuit attachment with or without inner cannula in place.



### ORDERING INFORMATION

Box quantity: 1

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

## PRODUCT BENEFITS

- Tubes are available in five color-coded adult sizes. Colour-coding throughout the D.I.C.® line enables caregivers to quickly identify and size tubes and related inner cannula

## PRODUCT SPECIFICATIONS

### Replacement Inner Cannula

Reorder Code	Tube ID (mm)	units / case
--------------	--------------	--------------

526060	5.0mm	60
--------	-------	----

526070	6.0mm	60
--------	-------	----

526080	7.0mm	60
--------	-------	----

526090	8.0mm	60
--------	-------	----

526100	9.0mm	60
--------	-------	----

### Uncuffed Flex D.I.C.® Tracheostomy Tubes

Reorder Code	Size	ID (mm)	OD (mm)	Length (mm)	Colour Code
--------------	------	---------	---------	-------------	-------------

504060	6.0 mm	6.0 mm	8.2 mm	64.0 mm	Orange
--------	--------	--------	--------	---------	--------

504070 7.0 mm 7.0 mm 9.6 mm 70.0 mm Green  
504080 8.0 mm 8.0 mm 10.9 mm 73.0 mm White  
504090 9.0 mm 9.0 mm 12.3 mm 79.0 mm Blue  
504100 10.0 mm 10.0 mm 13.7 mm 79.0 mm Yellow

Notes: Flex D.I.C.® tracheostomy tubes come with an obturator, Velcro® brand tube holder, and one colour-coded inner cannula.