

MOBILITY SCOOTER CATALOG



Get **Out**
There!

PRIDE® **MOBILITY SCOOTERS**
LIVE YOUR BEST®

pridemobility.com





Don't let a change in mobility affect your lifestyle. Live your best with Pride® Mobility Scooters and enjoy sleek, sporty style in a high-performance package.

With an array of models and colors to choose from, Pride Mobility Scooters make it easier for you to get out there and enjoy life.

What's Inside

- Go-Go® Travel Mobility 3-5
- Jazzy® Zero Turn 6
- Revo™ 2.0 6
- Victory® Series 7-8
- Pursuit® Series 9
- Maxima 10
- Raptor 10
- Scooter Accessories 11
- Specifications 13-14



Choosing the right Pride® Mobility Scooter

When choosing a Pride® Mobility Scooter, specific requirements and conditions can affect your process.

Indoor maneuverability

- Travel mobility to mid-size luxury scooters offer enhanced maneuverability

Stability

- 4-wheel scooters provide increased stability

Indoor and outdoor versatility

- Mid-size and full-size luxury scooters are great for both indoor and outdoor use

Outdoor capabilities

- Mid-size and full-size luxury scooters offer excellent outdoor capabilities

Increased weight capacity

- Mid-size and full-size luxury scooters provide increased weight capacity

Easy portability

- Travel mobility scooters are easy to disassemble, transport and store

Speed and performance

- Mid-size and full-size luxury scooters offer increased speed and performance

pridemobility.com



Pride®
Mobility Products Corp.
— LIVE YOUR BEST —



GO CHAIR[®]

Features

- 300 lbs. weight capacity
- One hand feather-touch disassembly
- 18" x 17", 5-point swivel seat with latch release
- Faster speed - up to 3.7 mph
- Compact turning radius 25.5"
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- Heaviest piece is only 36 lbs.



Feather-touch disassembly

Color Choices



Ruby Red

Sapphire Blue

Rose Quartz

Citrine Yellow

Amber Orange

Pearl White

3-WHEEL SCOOTER

GO GO[®] es²

Features

- 250 lbs. weight capacity
- Two speed settings - high and low (indoor/outdoor), maximum speed up to 4 mph
- Up to 6.8 miles per charge range
- Lightweight, one-piece frame for added stability and easy portability

Color Choices



Black



3-wheel



4-wheel

3 & 4-WHEEL SCOOTER

GO GO[®] ULTRAX

Features

- 260 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4 mph
- Up to 6.9 miles per charge range (3W), 7.2 per charge range (4W)

Color Choices



Red

3 & 4-WHEEL SCOOTER

GOGO
ELITE TRAVELLER®

Features

- 300 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4 mph
- Up to 10 miles per charge range (3W), 9.7 per charge range (4W)
- Charger port in tiller

Color Choices

INTERCHANGEABLE SHROUDS⁷



Red



Blue



3-wheel



4-wheel

3 AND 4-WHEEL SCOOTER

GOGO
ELITE TRAVELLER® PLUS

Features

- 300 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4.5 mph
- Up to 14 miles per charge range (3W), 8.3 per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS⁷



Red



Blue



3-wheel



4-wheel

3 AND 4-WHEEL SCOOTER

GOGO
SPORT

Features

- 325 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4.7 mph
- Up to 13 miles per charge range (3W), 12 miles per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS⁷



Red



Blue



3-wheel



4-wheel



3-wheel



4-wheel

3 AND 4-WHEEL SCOOTER

GOGO[®]


LX WITH CTS SUSPENSION

Features

- 300 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 4.4 mph
- Up to 9.4 miles per charge range (3W), 11 miles per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS⁷



3-WHEEL SCOOTER

GOGO[®]

FOLDING SCOOTER

Features

- 250 lbs. weight capacity
- Maximum speed up to 4 mph
- Up to 13 miles per charge range (lithium ion), 9.3 miles per charge range (sealed lead acid)
- Height-adjustable tiller
- Lighting package
- Optional lithium ion battery
- Optional armrests

Color Choices



White



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go[®] Folding Scooter and its lithium-ion battery is permitted.

NEW

4-WHEEL SCOOTER

Jazzy® ZERO TURN



Patented, intelligent-turning



Features

- 300 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- iTurn Technology™ provides 4-wheel stability with 3-wheel maneuverability
- 38" turning radius navigates tight corners and small spaces effortlessly
- Dual motors for added traction in rough terrain
- Standard under seat storage
- Delta tiller with wraparound handles
- Charger port in tiller
- LED lighting package
- USB charger built into the tiller for convenient charging of smart phone and portable devices

Color Choices

INTERCHANGEABLE SHROUDS*



NEW

3 AND 4-WHEEL SCOOTER

REVO™ 2.0



3-wheel




4-wheel

Features

- 375 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- CTS Suspension (Comfort-Trac front and rear suspension)
- Up to 12 miles per charge range (3W), 13 miles per charge range (4W)
- Standard under seat storage
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight
- Easily charge the batteries on or off the unit
- USB mobile device charging port

Color Choices

INTERCHANGEABLE SHROUDS*



VICTORY[®] SERIES



3-wheel



4-wheel

3 AND 4-WHEEL SCOOTER

VICTORY[®] 9

Features

- 300 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.3 mph
- Up to 13 miles per charge range
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port

Color Choices



Candy
Apple Red



Viper
Blue



3-wheel



4-wheel

3 AND 4-WHEEL SCOOTER

VICTORY[®] 10

Features

- 400 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.3 mph
- Up to 16 miles per charge range
- Delta tiller with wraparound handles
- Long lasting, high intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port

Color Choices



Candy
Apple Red



Viper
Blue

The Ultimate in Style and Performance®

4-WHEEL SCOOTER

VICTORY® 10 LX WITH CTS SUSPENSION

Features



- 400 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 5.6 mph
- Up to 16 miles per charge range
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- Full LED lighting package that includes curb and rear lights and directional signals
- Easy access tiller-mounted charger port

Color Choices



Candy
Apple Red



True
Blue

4-WHEEL SCOOTER

VICTORY® SPORT

Features



- 350 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 8 mph
- Up to 11 miles per charge range
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination with directional signals
- Front and rear suspension
- Easy access tiller-mounted charger port

Color Choices



Candy
Apple Red



Viper
Blue



FULL-SIZE LUXURY

*Outstanding Outdoor
Performance with
Luxurious Features*

4-WHEEL SCOOTER



Features

- 400 lbs. weight capacity
- Infinitely-adjustable tiller angle with conveniently located handle
- Maximum speed up to 8 mph
- Front and rear suspension
- Up to 14 miles per charge range (NF-22), 19 miles per charge range (G24)
- Feather-touch disassembly
- Delta tiller with wraparound handle
- Full lighting package that includes directional signals and lower headlight
- Reclining high-back seat with sliders and headrest

Color Choices



Candy Apple Red



Silver

4-WHEEL SCOOTER



Features

- 400 lbs. weight capacity
- Infinitely-adjustable tiller angle with conveniently located handle
- Maximum speed up to 8.5 mph
- Up to 30 miles per charge range (G27), 22 miles per charge range (G24)
- Front and rear suspension
- Delta tiller with wraparound handle
- Deluxe reclining high-back seat with sliders and headrest
- 13" low-profile, pneumatic tires for excellent outdoor performance

Color Choices



Candy Apple Red



Silver

RAPTOR

3-WHEEL SCOOTER

FULL-SIZE LUXURY

Outstanding Outdoor Performance with Luxurious Features



Features

- 400 lbs. weight capacity
- Front and rear suspension
- Maximum speed up to 14 mph
- Full LED lighting package with directional signals and brake lights that offer full visibility for increased safety
- Hand brake with regenerative braking system
- Secure, lockable under seat storage compartment with seat support rod
- Easy to read LCD console which shows estimated speed, mileage and battery display
- Limited recline fold-flat seat with manual adjustment and angle/width adjustable armrest
- Sleek and solid design for enhanced stability
- Front storage compartment located on the tiller
- Tubeless pneumatic tires for a smooth ride

Color Choices



Silver

MAXIMA

3 AND 4-WHEEL SCOOTER

HEAVY-DUTY

Heavy-Duty Performance at an Exceptional Value



3-wheel

Features

- 500 lbs. weight capacity
- Ultra heavy-duty drivetrain with 110 amp controller
- Maximum speed up to 5.3 mph
- Easily adjustable tiller
- Up to 18 miles per charge range (3W, NF-22), 9.6 miles per charge range (Power elevating seat, U-1), 16 miles per charge range (4W)
- Durable and stylish rear polyurethane bumper
- User-friendly console



4-wheel

Color Choices



Candy Apple Red



Viper Blue

ACCESSORIES

A wide selection of Pride Mobility Scooter accessories are available to personalize your Pride Mobility Scooter to better meet your needs and make it uniquely yours.



Cane clip attachment



Cup holder



Single cane holder



Safety flag



Double cane holder



Under seat storage bins



Rear basket



XLR USB charger



Oxygen holder



RAM® X-Grip®
cell phone holder



Wishbone crutch holder



Additional battery box
(12 amp)



Lapbelt



Saddlebag



Walker holder



Additional battery box
(18 amp)



Rearview mirror



Weather breaker



Weather cover

Military Patches Available



DISCLAIMER: Not all accessories are available for all products. Please speak with your Provider for details and availability.

INNOVATIVE TECHNOLOGY

NEW!



Patented, intelligent-turning iTurn Technology™ offers a **38" turning radius** to ensure effortless navigation within tight corners and small spaces.



CTS SUSPENSION



Front Suspension



Rear Suspension

Comfort-Trac Suspension (CTS) is the next generation in mobility scooter ride comfort. Pride's **industry-first** technology incorporates advanced suspension components for improved absorption over rough surfaces, and a smoother ride experience over various terrains.

FEATHER-TOUCH DISASSEMBLY



Scooter/travel mobility disassembly



Go Chair® disassembly

Most Go-Go® Travel Mobility and Pride® Mid-Size Scooters feature **Feather-Touch** — one-hand disassembly for hassle-free transport and storage.

SPECIFICATIONS

	Front Tires	Rear Tires	Maximum Speed ²	Ground Clearance ^{4,11,14}	3-W Turning Radius ⁴	4-W Turning Radius ⁴	3-W Length/Width ^{4,10}	4-W Length/Width ^{4,10}	Seat-to-Ground Height Range	(3/4-wheel) Weight ^{3,4,12}	
TRAVEL MOBILITY	Jazzy® Zero Turn	2"x7" SOLID	2.5"x9" SOLID	UP TO 4 MPH	1.3" at motor	N/A	38.25"	N/A	45.75"/20"	N/A	120 LBS.
	Revo™ 2.0	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	3" x 9" SOLID	UP TO 5 MPH	1.25" at transaxle	40"	52"	46"/21"	47"/21"	21-23.5"	120 LBS./131 LBS.
	Go Chair®	6" SOLID	9" SOLID	UP TO 3.7 MPH	2"	N/A	25.5" with foot platform	N/A	28.5" without footplate 33"/22.5"	22"-23"	67 LBS.
	Go-Go® ES 2 K0800 - GROUP 1 STANDARD	2" x 7.5" SOLID	3" x 7.5" SOLID	UP TO 4 MPH	1.5" at motor	31"	N/A	36.75"/19.5"	N/A	19"-20"	58 LBS.
	Go-Go® Ultra X	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4 MPH	1.25" at motor	33"	44"	37"/19.5"	39.5"/19.5"	18.25"-20.25"	76 LBS./80 LBS.
	Go-Go Elite Traveller® K0800 - GROUP 1 STANDARD	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4 MPH	2"	33"	44"	37"/19.5"	39.5"/19.5"	18.25"-20.25"	76 LBS./80 LBS.
	Go-Go Elite Traveller® Plus K0806 - GROUP 2 STANDARD ¹³	2" x 8" SOLID (3W)	3" x 9" SOLID	UP TO 4.5 MPH	1.5" at motor	34.5"	45.75"	40"/21.25"	40.25"/21.25"	20.5"-22"	86 LBS./87 LBS.
	Go-Go® Sport K0801 - GROUP 1 HEAVY DUTY	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	3" x 9" SOLID	UP TO 4.7 MPH	1.5" at motor	36"	50.5"	42"/21.25"	43"/21.25"	20.5"-22"	93 LBS./99 LBS.
	Go-Go® LX with CTS Suspension	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4.4 MPH	1.3" at motor	33"	44.75"	39.75"/20.5"	41"/20.5"	20.75"-22.25"	85LBS./87 LBS.
Go-Go® Folding Scooter	1.5" x 7" SOLID	2" x 8" SOLID	UP TO 4 MPH	2.3"	32"	N/A	37"/19"	N/A	20.5"	51 LBS.	
MID-SIZE LUXURY	Victory® 9 K0800 - GROUP 1 STANDARD	3W: 3" x 8" SOLID 4W: 3" x 9" SOLID	3" x 9" SOLID	UP TO 5.3 MPH	3W: 2" 4W: 2.25"	40.25"	51.5"	43.5"/22.25"	45.5"/22.5"	19.25"-20.75"	117 LBS./131 LBS.
	Victory® 10 K0801 - GROUP 1 HEAVY DUTY	3.5" x 10" SOLID	3.5" x 10" SOLID	UP TO 5.3 MPH	3" at mid-frame	42"	54"	45.25"/22.25"	47"/22.25"	20.5"-22"	133 LBS./137 LBS.
	Victory® 10 LX with CTS Suspension¹⁵	3.5" x 10" SOLID	3.5" x 10" SOLID	UP TO 5.6 MPH	2.25" at motor	N/A	60"	N/A	47"/22.25"	23.5"-25.5"	179 LBS.
	Victory® Sport¹⁵	3.5" x 10" SOLID	4" x 10" SOLID	UP TO 8 MPH	2.25"	N/A	51.75"	N/A	47"/26"	21.5"-24.25"	173 LBS.
HEAVY-DUTY	Maxima K0802 - GROUP 1 VERY HEAVY DUTY (3W) K0808 - GROUP 2 VERY HEAVY DUTY (4W)	3" X 10" SOLID	4" x 10" SOLID	UP TO 5.3 MPH	4.25" 3W 3.25" 4W	50.5"	54"	46"/25.5"	47"/25.5"	22.5"-24.5"	171 LBS./176 LBS.
FULL-SIZE LUXURY	Pursuit®¹⁵	3" X 13" SOLID	13" SOLID	UP TO 8 MPH	3.5"	N/A	67.75" LEFT 66.75" RIGHT	N/A	53.25"/25.25"	24.75"-26.25"	225 LBS.
	Pursuit® XL¹⁵	13" PNEUMATIC	13" PNEUMATIC	UP TO 8.5 MPH	5.25"	N/A	82.5"	N/A	56"/27.25"	26.25"-27.75"	254 LBS.
	Raptor¹⁵	3" X 16" PNEUMATIC	3" X 16" PNEUMATIC	UP TO 14 MPH	3.5" at center of frame	60"	N/A	59"/24.75"	N/A	28"	208 LBS.

- All shroud panels included.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Battery weight not included.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

- Battery pack is 18 lbs. with 12AH batteries.
- Deep cycle (AGM or gel cell type recommended).
- 18 AH battery pack weighs 29.5 lbs.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

Get Out There!

Seat Type (Standard)	Seat-to-Deck Height Range	Battery Weight ⁹	Battery Requirements ⁶	Per Charge Range (Max.) ^{2,8}	Battery Charger	Maximum Weight Capacity	Suspension	Weight of Heaviest Piece ⁴
FOLDABLE, MOLDED PLASTIC (17" X 17")	17.5"	13.5 LBS. EA.	(2) 22AH	UP TO 10 MILES	OFF-BOARD, 3.5A	300 LBS.	FRONT AND REAR	48 LBS. (REAR SECTION)
FOLDABLE, MOLDED PLASTIC (20" X 18")	18"	24 LBS. EA.	(2) U1	3W: UP TO 12 MILES 4W: UP TO 13 MILES	OFF-BOARD, 3.5A	375 LBS.	FRONT AND REAR	48 LBS. (3-WHEEL) (REAR SECTION) 57 LBS. (4-WHEEL) (FRONT SECTION)
COMPACT, FOLDABLE (18" X 17")	N/A	14 LBS. EA.	(2)18AH	UP TO 8.7 MILES	OFF-BOARD, 2A	300 LBS.	LIMITED	36 LBS. (REAR SECTION)
COMPACT, FOLDABLE (17" X 15")	16.25"-17.25"	9 LBS. EA.	(2)12AH	UP TO 6.8 MILES	OFF-BOARD, 1.5A	250 LBS.	N/A	41 LBS. (3-WHEEL) (SCOOTER BASE)
COMPACT, FOLDABLE (17" X 17")	15.75"-17.25"	9 LBS. EA. ⁵	(2)12AH	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	OFF-BOARD, 2A	260 LBS.	N/A	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)
COMPACT, FOLDABLE (17" X 17")	15.75"-17.25"	9 LBS. EA. ⁵ / 14 LBS. EA. ⁷	(2)12AH OR (2)18AH	3W: UP TO 7.2 MILES (12AH) ⁹ UP TO 10 MILES (18AH) ⁹ 4W: UP TO 6 MILES (12AH) ⁹ UP TO 9.7 MILES (18AH) ⁹	OFF-BOARD, 2A	300 LBS.	N/A	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)
COMPACT, FOLDABLE (17" X 17")	16.25"-17.75"	14 LBS. EA. ⁷	(2)18AH	3W: UP TO 14 MILES ⁸ 4W: UP TO 8.3 MILES ⁸	OFF-BOARD, 2A	300 LBS.	N/A	32 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)
COMPACT, FOLDABLE (17" X 17")	17"-18"	14 LBS. EA. ⁷	(2)18AH	3W: UP TO 13 MILES 4W: UP TO 12 MILES	OFF-BOARD, 2A	325 LBS.	N/A	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)
COMPACT, FOLDABLE (17" X 17")	16.25"-17.25"	9 LBS. EA. ⁵ / 14 LBS. EA. ⁷	(2)12AH OR (2)18AH	3W: UP TO 9.9 MILES (12AH) UP TO 9.44 MILES (18AH) 4W: UP TO 5.5 MILES (12AH) UP TO 11 MILES (18AH)	OFF-BOARD, 2A	300 LBS.	FRONT AND REAR	34 LBS. (3-WHEEL) (FRONT SECTION) 37 LBS. (4-WHEEL) (FRONT SECTION)
COMPACT, FOLDABLE (16" X 13")	13"	9 LBS. EA. (SEALED LEAD ACID) / 5 LBS. EA. (LITHIUM-ION)	(2)12AH OR (1) 12AH (LITHIUM-ION)	UP TO 13 MILES (LITHIUM-ION) UP TO 9.3 MILES (SEALED LEAD ACID)	OFF-BOARD, 2A (SEALED LEAD ACID) OFF-BOARD, 2.5A (LITHIUM-ION)	250 LBS.	N/A	51 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)
FOLDABLE, MOLDED PLASTIC (18" X 17")	17.25"-18.75"	24 LBS. EA	(2) U1	UP TO 13 MILES	OFF-BOARD, 3.5A (U-1)	300 LBS.	N/A	48 LBS. (3-WHEEL) (FRONT SECTION) 59 LBS. (4-WHEEL) (FRONT SECTION)
FOLDABLE, MOLDED PLASTIC (18" X 17")	17.5"-19"	24 LBS. EA. / 32 LBS. EA.	(2) U1 OR (2) 40AH	UP TO 16 MILES	OFF-BOARD, 3.5A (U-1) / OFF-BOARD, 5A (40AH)	400 LBS.	N/A	54 LBS. (3-WHEEL) (FRONT SECTION) 61 LBS. (4-WHEEL) (FRONT SECTION)
RECLINING H/B W/SLIDERS AND HEADREST (18" X 18")	16.5"-18.5"	32 LBS. EA.	(2) 40AH	UP TO 16 MILES	OFF-BOARD, 5A (40AH)	400 LBS.	FRONT AND REAR	78 LBS. (4-WHEEL) (FRONT SECTION)
RECLINING H/B W/SLIDERS AND HEADREST (18" X 18")	16"-18.5"	32 LBS. EA.	(2) 40AH	UP TO 11 MILES	OFF-BOARD, 5A	350 LBS.	FRONT AND REAR	68 LBS. (4-WHEEL) (FRONT SECTION)
DELUXE CONTOURED MEDIUM-BACK (22" X 18")	3W: 16"-18" 4W: 16"-18"	38 LBS. EA. (NF-22 OR 55AH) 24 LBS. EA. (U-1)	(2) NF-22 (2) U-1 (POWER ELEVATING SEAT)	3W: UP TO 18 MILES, U-1 (POWER ELEVATING SEAT) UP TO 9.6 MILES 4W: UP TO 16 MILES	OFF-BOARD, 5A	500 LBS.	3-W REAR / 4-W FRONT AND REAR	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)
RECLINING H/B W/SLIDERS AND HEADREST (18" X 18")	17"-18.75"	38 LBS. EA. (NF-22 OR 55AH) 53 LBS. EA. (75AH)	(2) NF-22 (2) GROUP 24	UP TO 14 MILES (NF-22) UP TO 19 MI. (GROUP 24)	OFF-BOARD, 5A (NF-22) OFF-BOARD, 8A (GROUP 24)	400 LBS.	FRONT AND REAR	95 LBS. (FRONT SECTION)
HIGH-BACK PILLOW TOP W/SLIDERS AND HEADREST (19.5" X 17.5")	21"-22.5"	53 LBS. EA. (75AH) 60 LBS. EA. (100AH)	(2) GROUP 27 OR (2) GROUP 24	G27: UP TO 30 MILES G24: UP TO 22 MILES	OFF-BOARD, 8A	400 LBS.	FRONT AND REAR	205 LBS. (MAIN FRAME)
HIGH-BACK (20" X 18") (100" LIMITED RECLINE)	19"	53 LBS. EA.	(2) 75AH	19 MILES (UP TO 400 LBS.) 31 MILES (UP TO 220 LBS. AND AT LOW SPEED)	OFF-BOARD, 8A	400 LBS.	FRONT AND REAR	N/A (UNIT DOES NOT DISASSEMBLE)

9. Battery weight may vary based on manufacturer.

10. Width of base (not including seat).

11. Does not account for anti-tip clearance.

12. Includes standard seat weight.

13. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.

14. At midframe.

15. Light package.

To find a compatible lift for your vehicle/mobility product, use our new

LIFT FINDER

pridemobility.com/lift-finder



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

AUTHORIZED PRIDE® PROVIDER

PROVIDER INFO

The information contained herein is correct at the time of publication;
we reserve the right to alter specifications without prior notice.

SCGUIDE-CAT-01-30-18