

# Sponges



PRM4408Z



PRM21408CZ



PRM21424

## Cotton Woven Gauze Sponges

Ideal for wound debriding, prepping, packing, dressing and general wound care.

- » 100% cotton sponges constructed of fine mesh gauze
- » Low-lint manufacturing process assures debris-free packaging
- » Not made with natural rubber latex

Item No.	Description	Pkg
<b>Sterile 2's</b>		
PRM2208Z	2 x 2" (5.08 cm x 5.08 cm), 8-Ply	2/pk, 50 pk/bx
PRM3312Z	3 x 3" (7.62 cm x 7.62 cm), 12-Ply	2/pk, 40 pk/bx
PRM4408Z	4 x 4" (10.16 cm x 10.16 cm), 8-Ply	2/pk, 50 pk/bx
PRM4412H	4 x 4" (10.16 cm x 10.16 cm), 12-Ply	2/pk

## Non-Sterile Bulk

PRM21208CZ	2 x 2" (5.08 cm x 5.08 cm), 8-Ply	200/sleeve
PRM21212CZ	2 x 2" (5.08 cm x 5.08 cm), 12-Ply	200/sleeve
PRM21312CZ	3 x 3" (7.62 cm x 7.62 cm), 12-Ply	200/sleeve
PRM21408CZ	4 x 4" (10.16 cm x 10.16 cm), 8-Ply	200/sleeve
PRM21412CZ	4 x 4" (10.16 cm x 10.16 cm), 12-Ply	200/sleeve

## Cotton Gauze Sponges

- » Not made with natural rubber latex

Item No.	Description	Pkg
<b>USP Type VII, Sterile 2's</b>		
PRM21420	2 x 2" (5.08 cm x 5.08 cm), 8-Ply	50 pk/bx, 30 bx/cs, 1500 pks/cs
PRM21422	3 x 3" (7.62 cm x 7.62 cm), 12-Ply	25 pk/bx, 48 bx/cs, 1200 pks/cs
PRM21424	4 x 4" (10.16 cm x 10.16 cm), 12-Ply	25 pk/bx, 24 bx/cs, 600 pks/cs
PRM21448	4 x 4" (10.16 cm x 10.16 cm), 8-Ply	25 pk/bx, 24 bx/cs, 600 pks/cs

## USP Type VII, Sterile 10's, Soft Pouch

PRM21426	4 x 4" (10.16 cm x 10.16 cm), 12-Ply	10 sponges/pk, 128 pk/cs, 1280/cs
PRM21428	4 x 4" (10.16 cm x 10.16 cm), 16-Ply	10 sponges/pk, 128 pk/cs, 1280/cs

## USP Type VII X-Ray Detectable, Sterile, Soft Pouch

PRM21430LF	4 x 4" (10.16 cm x 10.16 cm), 16-Ply	10 sponges/pk, 135 pk/cs, 1350/cs
------------	--------------------------------------	-----------------------------------