

## Pro-Vent™ Mask Configurations

Configuration	Part #	NDC #	HCPCS #
Bubbles the Fish™ II (Ages 3+ years)	#044F7248	#83490-0440-68	A7015
PARI LC® Adult Mask	#044F7251	#83490-0440-78	A7015


## Features

- Front-load design directs aerosol toward the patient's mouth
- Comfortable aerosol delivery with soft, flexible masks
- Pediatric character design makes treatments fun for kids!



## The Difference is in the Delivery

■ Bubbles™ II with PARI LC® Plus  
 ■ Salter Mask with Misty Neb



### Accessories:

PARI LC® Sprint = 023F35  
 PARI LC® Plus = 022F81  
 PARI Pals™ Activity Book = 044D0029  
 PARI Pals™ Stickers = 044D0030

**KEY FINDINGS: A higher Inhaled Mass combined with lower facial and eye deposition gives patients increased drug delivery to the lungs with reduced exposure to the eyes and skin.**

<sup>1</sup>Facial and Ocular Deposition of Nebulized Budesonide. K. Harris, G. Smaldone. CHEST 2008; 133:482-488. Tables 2-4.  
<sup>2</sup>All registered marks and trademarks belong to their respective companies.

**PRO Vent™**  
 TECHNOLOGY



Fun, safe, optimized delivery of aerosol.

Only the PARI adult mask and Bubbles The Fish™ II pediatric mask are designed with Pro-Vent™ technology, proven to deliver more medication while reducing facial and eye deposition.

*From the leader in clinically proven aerosol delivery*