

# RIB SUPPORTS



## MODEL: 2459



### RIB BELT FOR MEN HCPCS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Wraparound elastic web material with adjustable contact closure fastening.
- 6" deep all around

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

#### MEASURE THE CHEST CIRCUMFERENCE AT THE BASE OF THE BREAST BONE

2459-UR	UR	29 – 38 in / 73.6-96.5 cm
2459-UL	UL	39 – 48 in / 99-121.9 cm

## MODEL: 2659



### RIB BELT FOR WOMEN HCPCS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Wraparound elastic web material with adjustable contact closure fastening.
- Cut low in front to accommodate the bust line
- 4" in front, 6" in back

INDICATIONS: Compression and support for rib fractures and contusions around the rib cage

#### MEASURE THE CHEST CIRCUMFERENCE AT THE BASE OF THE BREAST BONE

2659-UR	UR	29 – 38 in / 73.6-96.5 cm
2659-UL	UL	39 – 48 in / 99-121.9 cm

# ABDOMINAL AND HERNIA SUPPORTS



## MODEL: 2509



### THREE PANEL ELASTIC BINDER FOR MEN HCPCS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Hinged three-panel elastic molds around surface anatomy for a comfortable, custom fit
- Can also be used as a rib fracture support
- 9" deep all around

INDICATIONS: Strains and weaknesses in the abdominal area, loss of muscle tone, post-surgical application

#### MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2509-S	S	29-32 in / 73.6-81.3 cm
2509-M	M	33-36 in / 81.3-91.4 cm
2509-L	L	37-40 in / 94-101.6 cm
2509-XL	XL	41-44 in / 104.1-111.8 cm
2509-2XL	2XL	45-48 in / 114.3-121.9 cm

## MODEL: 2689



### THREE PANEL ELASTIC BINDER FOR WOMEN HCPCS: A4466 • MINIMUM

- Form fitting support relieves discomfort through circular compression
- Hinged three-panel elastic molds around surface anatomy for a comfortable, custom fit
- Contoured to fit the female form
- 9" deep all around

INDICATIONS: Strains and weaknesses in the abdominal area, loss of muscle tone, post-surgical application

#### MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2689-S	S	19-22 in / 48.3-55.9 cm
2689-M	M	23-26 in / 58.4-66 cm
2689-L	L	27-30 in / 68.6-76.2 cm
2689-XL	XL	31-34 in / 78.7-86.4 cm
2689-2XL	2XL	35-38 in / 88.9-96.5 cm

## MODEL: 2518



### MULTIPLE USE ELASTIC ABDOMINAL BINDER HCPCS: A4466 • MINIMUM

- Lightweight full length binder has a variety of applications, most commonly to protect post-surgical dressings
- Effective in supporting abdominal contents when there has been a loss of muscle tone
- Flexible plastic support stays prevent rolling
- 10" deep all around

INDICATIONS: Umbilical hernia, weak or pendulous abdomen, abdominal contusions and tenderness, minor rib fracture

#### MEASURE AROUND THE FULLEST PART OF THE ABDOMEN

2518-S	S	29-32 in / 73.6-81.3 cm
2518-M	M	33-36 in / 81.3-91.4 cm
2518-L	L	37-40 in / 94-101.6 cm
2518-XL	XL	41-44 in / 104.1-111.8 cm
2518-2XL	2XL	45-48 in / 114.3-121.9 cm

## MODEL: 2956



### LIGHTWEIGHT BEIGE HERNIA SUPPORT HCPCS: L8300 (SINGLE) •• MEDIUM L8310 (BILATERAL)

- Plush lightweight material is contact closure compatible
- Elastic waistband in back molds around surface anatomy for a comfortable, custom fit
- Removable foam pads provide gentle, localized support at the site of the hernia
- Soft elastic understraps can be tightened or loosened to increase or lessen support
- Completely free of any metal component

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

#### MEASURE AROUND THE WIDEST PART OF THE HIPS

2956-S	S	30-35 in / 76-89 cm
2956-M	M	35-40 in / 89-102 cm
2956-L	L	41-46 in / 104-117 cm





## HERNIA SUPPORTS

MODEL: 2936  
SINGLE  
MODEL: 2939  
BILATERAL



**HERNIA SUPPORT (SCROTAL ELASTIC TRUSS)**  
HCPCS: L8300 (2936) •• MEDIUM  
L8310 (2939)

- Relieves discomfort while treating inguinal or scrotal hernias
- Design of the truss allows the specially shaped pad to self-align and maintain its proper position, thus applying support as a constant, unified force.

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

MEASURE AROUND THE  
WIDEST PART OF THE HIPS

2936/L-S	S	LEFT	30-32 in / 76.2-81.3 cm
2936/L-M	M	LEFT	34-36 in / 86.4-91.4 cm
2936/L-L	L	LEFT	38-40 in / 96.5-101.6 cm
2936/L-XL	XL	LEFT	42-44 in / 106.7-111.8 cm
2936/L-2XL	2 XL	LEFT	45-48 in / 114.3-121.9 cm
2936/R-S	S	RIGHT	30-32 in / 76.2-81.3 cm
2936/R-M	M	RIGHT	34-36 in / 86.4-91.4 cm
2936/R-L	L	RIGHT	38-40 in / 96.5-101.6 cm
2936/R-XL	XL	RIGHT	42-44 in / 106.7-111.8 cm
2936/R-2XL	2 XL	RIGHT	45-48 in / 114.3-121.9 cm
2939-S	S	BILATERAL	30-32 in / 76.2-81.3 cm
2939-M	M	BILATERAL	34-36 in / 86.4-91.4 cm
2939-L	L	BILATERAL	38-40 in / 96.5-101.6 cm
2939-XL	XL	BILATERAL	42-44 in / 106.7-111.8 cm
2939-2XL	2 XL	BILATERAL	45-48 in / 114.3-121.9 cm

MODEL: 2958



**LIGHTWEIGHT SINGLE HERNIA SUPPORT**  
HCPCS: L8300 •• MEDIUM

- Holds a single inguinal or scrotal hernia comfortably in reduction
- Can also be used in special situations after a hernia repair to prevent recurrence
- Low profile asymmetric design allows it to be worn under clothing without detection
- 4" web elastic encircling belt molds around body contours for a comfortable, custom fit

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

MEASURE AROUND THE WIDEST  
PART OF THE HIPS

2958-S	S	30-32 in / 76.2-81.3 cm
2958-M	M	34-36 in / 86.4-91.4 cm
2958-L	L	38-40 in / 96.5-101.6 cm
2958-XL	XL	42-44 in / 106.7-111.8 cm

Specify for hernia on wearer's **RIGHT** or **LEFT** side

MODEL: 2958D



**LIGHTWEIGHT DOUBLE HERNIA SUPPORT**  
HCPCS: L8310 •• MEDIUM

- Relieves discomfort while treating inguinal or scrotal hernias
- Can also be used after hernia repair to prevent recurrence
- Removable foam pads provide localized support at hernia site
- Detachable understraps can be tightened or loosened to increase or lessen support

INDICATIONS: Reducible inguinal / scrotal hernia, single or double

MEASURE AROUND THE WIDEST  
PART OF THE HIPS

2958D-S	S	30-32 in / 76.2-81.3 cm
2958D-M	M	34-36 in / 86.4-91.4 cm
2958D-L	L	38-40 in / 96.5-101.6 cm
2958D-XL	XL	42-44 in / 106.7-111.8 cm



## CERVICAL AND CLAVICLE SUPPORTS

MODEL: 2452



**POSTURE SUPPORT**  
HCPCS: A4466 •• MINIMUM

- Serves as a gentle reminder to hold the shoulders in a natural, correct posture, especially while performing highly repetitive tasks like data entry
- Also helps prevent strain and fatigue
- Made from finest quality beige hinged elastic web: 3" wide shoulder straps, 6" back panel

INDICATIONS: Slumping shoulders, poor posture

MEASURE WAIST CIRCUMFERENCE

2452-XS	XS	25-28.5 in / 63-72 cm
2452-S	S	28.5-31.5 in / 72-80 cm
2452-M	M	31.5-35.5 in / 80-90 cm
2452-L	L	35.5-39.5 in / 90-100 cm
2452-XL	XL	39.5-43.5 in / 100-110 cm

MODEL: 2455



**POSTURE CONTROL SHOULDER BRACE**  
HCPCS: L3650 •• MEDIUM

- Serves as a gentle reminder to hold the shoulders in a natural, correct posture, especially while performing highly repetitive tasks like data entry
- Also helps prevent strain and fatigue
- Comfortable white canvas back panel, non-elastic web shoulder straps

INDICATIONS: Slumping shoulders, poor posture

MEASURE AROUND CHEST AT THE  
BASE OF THE BREAST BONE

2455-XS	XS	UP TO 26 in / UP TO 66 cm
2455-S	S	27-32 in / 68.9-81.3 cm
2455-M	M	33-38 in / 83.8-96.5 cm
2455-L	L	39-44 in / 99-111.7 cm
2455-XL	XL	OVER 45 in / OVER 114.3 cm

