

One Family Fits All

Introducing the Heelift® Family

Now there is a Heelift® Suspension Boot to help *all* of your at-risk patients prevent costly and painful heel pressure ulcers. With the heel being the second most common site where pressure ulcers develop, you cannot be too careful. And now that CMS has stated that it will no longer reimburse for hospital-acquired pressure ulcers, prevention is a must.

DM Systems has answered your requests for a smaller Heelift for your Petite/Pediatric patients and a much larger model for your Bariatric patients. The two new sizes share all of the same features that have made Heelift Suspension Boot the proven choice for hundreds of hospitals, clinics and long term care facilities.



Standard*

Bariatric

Petite*

There is only one Heelift®.

*Available in either convoluted or smooth foam.

Heelift® Standard

Available in Smooth and Convoluted Foam

Our signature model has been updated and enhanced numerous times over the years based on our research and your recommendations. Includes the extra pad which can be used to help prevent hip rotation or foot drop. The pad can easily be customized to address Achilles tendon erythema.



Indicated for Patients:

Calf Circumference	Height Range	Weight Range
8 inches – 14 inches 20.32 cm – 35.56 cm	60 inches – 77 inches 152.40 cm – 195.58 cm	120 lbs. – 250 lbs. 54.43 kilo – 113.40 kilo

Heelift® Petite

Available in Smooth and Convoluted Foam

All of the same features of the Heelift® Standard but in the perfect size to accommodate smaller adults or pediatric patients.



Indicated for Patients:

Calf Circumference	Height Range	Weight Range
6 inches – 10 inches 15.24 cm – 25.40 cm	46 inches – 64 inches 116.84 cm – 162.56 cm	70 lbs. – 130 lbs. 31.75 kilo – 58.97 kilo

Heelift® Bariatric

Available in Smooth Foam

All of the same features of the Heelift® Standard but in a larger size and with additional thickness to address the needs of bariatric patients.



Indicated for Patients:

Calf Circumference	Height Range	Weight Range
12 inches – 23 inches 30.48 cm – 58.42 cm	65 inches – 80 inches 165.10 cm – 203.20 cm	220 lbs. – 600 lbs. 99.79 kilo – 272.16 kilo

All Heelift® sizes feature:



- Extra pad which can be used to help prevent foot drop or hip rotation



- The pad can easily be customized to address Achilles tendon erythema



- Friction-free tricot backing for improved patient mobility, and a Polyurethane stiffener to prevent buckling



- Non-abrasive, easy-open straps with D-ring closures allow easy, one-handed adjustment of straps for ultimate comfort



- Extended stitching to narrow the forefoot and increase support
- Latex-free
- Machine washable



1316 Sherman Avenue
Evanston, IL 60201 USA
800.254.5438 • Fax: 847.328.9561
www.heelift.com
info@dmsystems.com

Heelift Standard, Heelift Petite and Heelift Bariatric are all covered by Patent No. 5449339 Suggested Code: E0191