


Instant Cold Packs

Cardinal Health™ Instant Cold Packs help provide consistent cold therapy for convenient and effective pain relief. Our comprehensive portfolio includes a variety of different options and sizes, including standard, insulated and premium non-sweat products. They are made with materials that do not contain DEHP.


#1
Instant Cold Pack
in hospitals!




Activation
is as easy as 1, 2, 3



1
Squeeze sides firmly where indicated to activate



2
Shake pack well to mix contents



3
Apply to desired area for 30 minutes

Features and benefits:

- Consistent and effective cold therapy
- Easy to activate — just squeeze and shake
- Can be applied directly to skin — no additional cover needed

Ordering information

Cat. No.	Description	Size (in.)	Packaging
103	Instant Cold Pack, Extra Small	5 x 5½ in.	50 ea/cs
103B	Instant Cold Pack, Extra Small, Bulk	5 x 5½ in.	50 ea/cs
102	Instant Cold Pack, Medium	6 x 6½ in.	16 ea/cs
11440-900	Instant Cold Pack, Large	6 x 9 in.	16 ea/cs
20204	Instant Cold Pack, Insulated, Medium	6 x 6½ in.	24 ea/cs
20104	Instant Cold Pack, Insulated, Large	6 x 9 in.	24 ea/cs
11440-512	Instant Cold Pack, Non-Sweat, Small	4½ x 9 in.	12 ea/bx, 4 bx/cs, 48 ea/cs
11445-020	Instant Cold Pack, Non-Sweat, Medium	6 x 6½ in.	20 ea/bx, 4 bx/cs, 80 ea/cs
11440-012	Instant Cold Pack, Non-Sweat, Large	6 x 9 in.	12 ea/ bx, 4 bx/cs, 48 ea/cs