

IV
Solutions

Enteral
Care

Respiratory
Care

Medical/
Surgical

Urology



The **Dale Nasal Dressing Holder** replaces adhesive tape to hold nasal dressings in place following rhinoplasty, septoplasty, sinus surgery, facial trauma, and nose bleeds.

Dale[®]
Always Better

Dale
600
CE +1-508-695-9316 REV G



Nasal Dressing Holder

REPLACES ADHESIVE TAPE TO HOLD NASAL DRESSINGS IN PLACE FOLLOWING RHINOPLASTY, SEPTOPLASTY, SINUS SURGERY, FACIAL TRAUMA AND NOSE BLEEDS.

It consists of a clear plastic Dressing Holder, and is secured by ear loops. Its hook and loop fastener tabs can be infinitely adjusted to fit any patient comfortably. The soft material of the Dressing Holder adds to the comfort, while the plastic center section holds dressings in place without absorbing secretions. The Nasal Dressing Holder is easy to remove and replace, making dressing changes simple and painless. And because it requires no adhesive tape, the Nasal Dressing Holder won't damage skin.

- Eliminates irritation and excoriation of facial skin
- Simplifies dressing changes
- Saves nursing time
- Improves patient comfort
- One size fits most patients
- Not made with natural rubber latex.

ORDERING INFORMATION

Product Number	Size	Quantity
600	One size fits most	10 per box

Dale[®]
www.dalemed.com

Dale[®] is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
©2017 Dale Medical Products, Inc.

PS-600 REV D