



## PROCARE®

### Quick-Fit® Basic Knee Splint

Foam laminate construction with 3" elastic contact closure straps provides compressive support and immobilization. Dual contoured, I-beam posterior stays and moveable medial/lateral double stay casings for correct fit.

**Indication:** Ideal for immobilization of the knee joint.

	Length	Max Thigh Circum.*	Size
79-96011	10"	29"	Univ.
79-96012	12"	29"	Univ.
79-96014	14"	29"	Univ.
79-96016	16"	29"	Univ.
79-96018	18"	36"	Univ.
79-96019	20"	36"	Univ.
79-96021	22"	36"	Univ.
79-96024	24"	36"	Univ.
79-96026	26"	36"	Univ.

\*Circumference 4" above patella.