



## Urine Collection Bags & Bottles

### Urine Collection Bag & Bottle Replacement Components/Parts

#### Urocare® Urinary Drainage Bottle Cap


Large, white cap. Ribbed for a sure grip and easy opening. Replacement cap for Urocare® Urinary Drainage Bottle on page 5. Anti-Reflux valve sold separately—see item #6001 below for more details.

Sold by the each, 1 cap per package.

This is a stocked item.

Non-Sterile, Latex-Free

**LATEX-FREE**



LEARN MORE



Item No.	Description	Dimensions†
<b>4101</b>	Bottle Cap	1.1" x 3.63" DIA. (2.8cm x 89mm Ø)

Case Item No.  
**410112** (12/CS)

#### Urocare® Urinary Drainage Bottle Hooks


Replacement hooks for Urocare® Urinary Drainage Bottle. Made from high-density polyethylene (HDPE), these hooks are strong yet light-weight and won't corrode or rust.

Sold by the pair, 1 pair per package.

This is a stocked item.

Non-Sterile, Latex-Free

**LATEX-FREE**



LEARN MORE



Item No.	Description	Dimensions†
<b>4102</b>	Bottle Hooks	3.63" x 1.9" x 0.22"

Case Item No.  
**410212** (12/CS)

#### Urocare® Gasket-Rings, Small

Replacement gasket-ring for: Urocare® Standard Bottom Drain Valve (#6004), Straight-Thru Adaptor (#6005), Male Urinal Adaptor (#6022), and all Quick Drain Valves™ (#700110, #700118 & #700172).

Sold by the package, 10 rings per package.

This is a stocked item.

Non-Sterile, Contains Latex



LEARN MORE



Item No.	Description	Dimensions†	Case Item No.
<b>5999</b>	Gasket-Rings, Small	3" x 5" x 0.25"	Not Available

Case Item No.  
Not Available

#### Urocare® Gasket-Rings, Large

Since Urocare's "High-Flow" Little Red Valve™ serves not only as an anti-reflux valve, but as a gasket-ring as well, in circumstances where the valve is not used, this gasket-ring must be used in place of the valve to prevent leaks in all Urocare® Reusable Latex Leg Bags and Urocare® Urinary Drainage Bottle which incorporate this valve.

Sold by the package, 10 rings per package.

This is a stocked item.

Non-Sterile, Contains Latex



LEARN MORE



Item No.	Description	Dimensions†	Case Item No.
<b>6000</b>	Gasket-Rings, Large	3" x 5" x 0.25"	Not Available

Case Item No.  
Not Available

#### Urocare® "High-Flow" Little Red Valve™

This new, **high-flow, self-cleaning**, patented, anti-reflux valve **prevents back-flow** while allowing fast, uninterrupted drainage for problem-free use. Don't settle for urine reflux or back-pressure caused by inferior flutter valves. Nothing beats Urocare's new "High-Flow" Little Red Valve™ for performance, reliability and safety. Replacement anti-reflux valve for use in all Urocare® Reusable Latex Leg Bags and Urocare® Urinary Drainage Bottle.

Sold by the each, 1 high-flow anti-reflux valve per package.

This is a stocked item.

Non-Sterile, 100% Silicone

**LATEX-FREE**



LEARN MORE



Item No.	Description	Dimensions†
<b>6001</b>	Little Red Valve™	0.65" DIA. x 0.77"

Case Item No.  
**600110** (10/CS)

† Approximate dimensions in inches(″) and centimeters(cm).