

Home Bed Assist Rail Folding Bed Board Combo




FIGURE A



FIGURE B



FIGURE C



RETAIL BOX

- Assistance for getting in and out of bed
- Can be used on either side of the bed (Figure A)
- Tri-fold bed board for easy portability and storage (Figure B)
- Removable handle
- Easy, quick, no tool assembly
- Comes in two pieces
- Simply slides under the mattress for extra-firm support (Figure C)
- Removable powder coated steel handle enables user to remove when not in use
- Handle with the "Mid Bar", creates grip for any height use
- Retail Packaged

PRODUCTS	Description
15062	1/cs

SPECIFICATIONS	Description
Rail Dimensions	11.5" (W) x 22.5" (H) x 30" (D)
Board Dimensions	26" (L) x 15" (W) x ¾" (H)
Weight Limit	250 lbs.
Warranty	Limited Lifetime