

COMPRESSION
THERAPY
COLLECTION



TRUFORM[®]
LEG HEALTH



Contents

| | |
|---------------------|---|
| Compression Profile | 3 |
| Therapy Guide | 4 |
| Measuring | 5 |
| Merchandising | 6 |

Fashion Wear

| | |
|------------------------|----|
| Ladies' Pattern Sheer | 8 |
| Ladies' Sheer | 10 |
| Ladies' Opaque | 14 |
| Ladies' Pattern Socks | 18 |
| Men's Dress Socks | 20 |
| Men's & Ladies' Casual | 22 |

Medical Leg Wear

| | |
|-----------------|----|
| Classic Medical | 24 |
| Anti-Embolism | 30 |

Specialty Wear

| | |
|-------------------------|----|
| Diabetic & Health Socks | 32 |
| Upper Extremity | 35 |

| | |
|---------------------|----|
| TRUFORM Accessories | 36 |
|---------------------|----|

| | |
|-------------------------|----|
| TRUFORM Product Summary | 38 |
|-------------------------|----|

Compression Profile

TRUFORM's 6 Points of Compression

To be therapeutically effective, medical stockings must provide a pressure reduction from the ankle upwards.

Truform takes compression measurements at precise intervals between the ankle and proximal end of the stocking - 4 points for below knee, 6 points for thigh length and pantyhose styles - to insure an accurate distribution of the pressure gradient.

Multiple measurements guarantee that Truform garments are comfortable and will deliver the prescribed graduated elastic compression therapy, without undesirable pressure inversions, when properly sized and fit.



Proximal Thigh ▶


Mid Thigh ▶

Fibular Head ▶

Mid Calf ▶

Lower Calf ▶

Ankle ▶



Therapy Guide

There is an appropriate Truform elastic compression therapy product for each of your customers.

Mild: 8-15 mmHg

Leg Fatigue
Slight ankle, leg and foot swelling

Moderate: 15-20 mmHg

Minor varicosities
Minor varicosities during pregnancy
Tired legs
Minor ankle, leg and foot swelling
Post surgery

Anti-Embolism: 18 mmHg

Prevention of thromboembolism
Reclining or bed confined home convalescence


Firm: 20-30 mmHg

Tired, aching legs
Moderate varicosities
Moderate varicosities during pregnancy
Moderate edema
Post sclerotherapy and other surgery
Help prevent recurrence of venous ulcerations
Superficial thrombophlebitis

Extra Firm: 30-40 mmHg

Severe varicosities
Severe edema
Lymphatic edema
Helps prevent and manage venous ulcerations
Chronic venous insufficiency
Helps prevent post thrombotic syndrome
Manage manifestations of PTS
Orthostatic hypotension
Post sclerotherapy
Post surgery

Compression Packaging Label



Available styles



Knee High



Thigh High




Pantyhose

All compression ranges stated are the mean compression for an average ankle for a specific size. Please refer to the product insert, product label and package for comprehensive instructions on the proper selection and safe use of Truform products. Sizing for Truform support hosiery is determined by combinations of circumference and length measurements, or by combinations of height and weight data. Some items, due to their design, may be "quick-fit" by shoe or sock size. Sizing charts for Truform hosiery are located on each product page throughout this catalog as appropriate.

MEASURING


Please follow these instructions carefully to determine the correct size for different product styles. Measurements should always be taken in the following order, and as soon as possible after waking when swelling is at its minimum.



1 Ankle circumference



2 Calf circumference



3 Knee length from the right-angle bend of the knee to the floor



4 Widest circumference of the thigh



5 Leg length from the floor up to the mid-thigh (1) or to the gluteal fold (2) for Anti-Embolism products.

Scan this QR Code to view the TRUFORM fitting video:

TRUFORM
LEG HEALTH

TRUFORM
LEG HEALTH


General Guidelines

- Measurements should always be taken as soon as possible after waking when swelling is at its minimum
- Measure against bare skin
- The measuring system is based on having the foot flat against the floor and leg at a right angle to the foot
- Measurements must be taken by pulling the tape snugly around the leg but without constricting the skin



TRUFORM
LEG HEALTH

Truform's planogram layouts make it easier for you to merchandise, for your fitters to sell from, and for your more self-service clients to shop. Modular planogram components enable you to grow your business gradually, purchasing at your own pace, and getting the most return for your investment. Choose between starter and intermediate assortments as well as product line specific assortments, depending on your business' particular needs and focus.





Truform products exhibit a dynamic contemporary package that educates the consumer and creates a bold retail presence. Coupled with POS (Point of Sale) displays, the Truform compression collection drives shelf impact to create distinctiveness and appeal. Innovative package design features create planograms that are easy to shop and easy to merchandise.



Ladies' Pattern Sheer



1782 15-20 mmHg
Knee High
BLACK



Ladies' Pattern Sheer Knee High

Closed Toe

- 1782 Dot Pattern 15-20 mmHg
- 1783 Diamond Pattern 15-20 mmHg



1783 BLACK



1782 BLACK

Material Information

1782 | 1783
 Manufactured using 80% nylon and 20% spandex

Product Description


Provide moderate graduated compression to help energize the legs, and reduce leg discomfort. Truform's Pattern Sheers provide long-lasting support that appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Moderate support
- Fashionable sheer fabric
- Available in diamond or dot patterns
- Knit-in heel and soft band for proper fit and comfort

LATEX FREE

Available Colors



| | CIRCUMFERENCE | | |
|---------|---------------|----------------|--------------|
| | ANKLE | CALF | |
| SMALL | 7" - 8.25" | 11" - 15" | 1782 1783 |
| MEDIUM | 8.25" - 9.5" | 11.75" - 16.5" | |
| LARGE | 9.5" - 10.75" | 12.5" - 18" | |
| X-LARGE | 10.75" - 12" | 13" - 19" | |

| Ladies' Pattern Sheer Knee High Compression Stockings | | | |
|---|--------------|--------------------------|------------|
| Model | Compression | Description | Latex Free |
| 1782 | 15 - 20 mmHg | Knee High Dot Pattern | ✓ |
| 1783 | 15 - 20 mmHg | Knee High Diamond Patten | ✓ |

Ladies' Sheer

1775 15-20 mmHg
Pantyhose BLACK



Ladies' Sheer Knee High

Closed Toe

- 1763 LITES 8-15 mmHg
- 1773 LITES 15-20 mmHg
- 0263 TruSheer 20-30 mmHg
- 0253 TruSheer 30-40 mmHg

Open Toe

- 1772 LITES 15-20

1763 BLACK

1772 NUDE

Material Information

1763
Manufactured using 88% nylon and 12% spandex

1772 | 1773
Manufactured using 80% nylon and 20% spandex

0263 | 0253
Manufactured using 57% nylon and 43% spandex

Product Description

Provide mild, moderate, firm, or extra firm graduated compression to help energize the legs, and reduce leg discomfort. They provide long-lasting support that appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm, or Extra-Firm support
- Fashionable sheer fabric
- Knit-in heel (1772, 1773, 0253, 0263)

Available Colors



*Colors may not be available in all styles.
Please refer to the Product Summary on pages 38 and 39 for more information.

LATEX FREE

| Ladies' Sheer Knee High Compression Stockings | | | |
|---|--------------|--------------------|------------|
| Model | Compression | Description | Latex Free |
| 1763 | 8 - 15 mmHg | Knee High | ✓ |
| 1772 | 15 - 20 mmHg | Knee High Open Toe | ✓ |
| 1773 | 15 - 20 mmHg | Knee High | ✓ |
| 0263 | 20 - 30 mmHg | Knee High | ✓ |
| 0253 | 30 - 40 mmHg | Knee High | ✓ |

| SHOE SIZE | | 1763 |
|-----------|---------------|------|
| SMALL | 4 - 5 | |
| MEDIUM | 5 1/2 - 7 1/2 | |
| LARGE | 8 - 10 1/2 | |
| X-LARGE | 11 - 12 1/2 | |

| CIRCUMFERENCE | | | 0263 0253 1772 1773 |
|---------------|---------------|----------------|------------------------------|
| ANKLE | | CALF | |
| SMALL | 7" - 8.25" | 11" - 15" | |
| MEDIUM | 8.25" - 9.5" | 11.75" - 16.5" | |
| LARGE | 9.5" - 10.75" | 12.5" - 18" | |
| X-LARGE | 10.75" - 12" | 13" - 19" | |

| CIRCUMFERENCE | | | 1773 |
|---------------|----------------|-------------|------|
| ANKLE | | CALF | |
| 2X-LARGE | 12.5" - 13.75" | 19.5" - 22" | |
| 3X-LARGE | 12.5" - 13.75" | 21.5" - 24" | |

Ladies' Sheer Thigh High

Closed Toe

● 1764 LITES 8-15 mmHg

● 1774 LITES 15-20 mmHg

● 0264 TruSheer 20-30 mmHg

● 0254 TruSheer 30-40 mmHg

Silicone stay-up top



Silicone stay-up top



264 BLACK

Material Information

1764
Manufactured using 88% nylon and 12% spandex

1774
Manufactured using 83% nylon and 17% spandex

0264 | 0254
Manufactured using 57% nylon and 43% spandex

Product Description


Provide mild, moderate, firm, or extra firm graduated compression to help energize the legs, and reduce leg discomfort. They provide long-lasting support that appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm, or Extra-Firm support
- Fashionable sheer fabric
- Attractive lace stay-up top

LATEX FREE

Available Colors



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

Ladies' Sheer Thigh High Compression Stockings

| Model | Compression | Description | Latex Free |
|-------|--------------|-------------|------------|
| 1764 | 8 - 15 mmHg | Thigh High | ✓ |
| 1774 | 15 - 20 mmHg | Thigh High | ✓ |
| 0264 | 20 - 30 mmHg | Thigh High | ✓ |
| 0254 | 30 - 40 mmHg | Thigh High | ✓ |

CIRCUMFERENCE

| | ANKLE | CALF | THIGH |
|---------|---------------|----------------|-----------------|
| SMALL | 7" - 8.25" | 11" - 15" | 15.75" - 24.5" |
| MEDIUM | 8.25" - 9.5" | 11.75" - 16.5" | 18" - 27.5" |
| LARGE | 9.5" - 10.75" | 12.5" - 18" | 21.25" - 30.75" |
| X-LARGE | 10.75" - 12" | 13" - 19" | 23.5" - 32" |

0264
0254

CIRCUMFERENCE

| | ANKLE | CALF | THIGH |
|---------|---------------|-------------|-------------|
| SMALL | 7" - 8.25" | 11" - 14" | up to 18.5" |
| MEDIUM | 8.25" - 9.75" | 13.5" - 16" | up to 19.5" |
| LARGE | 9.75" - 11" | 15.5" - 18" | up to 20.5" |
| X-LARGE | 11" - 12.5" | 17.5" - 20" | up to 21.5" |

1774

CIRCUMFERENCE

| | SHOE SIZE | THIGH |
|---------|---------------|-------|
| SMALL | 4 - 5 | 18.5" |
| MEDIUM | 5 1/2 - 7 1/2 | 19.5" |
| LARGE | 8 - 10 1/2 | 20.5" |
| X-LARGE | 11 - 12 1/2 | 21.5" |

1764

Ladies' Sheer Pantyhose

Closed Toe

- 1765 LITES 8-15 mmHg
- 1775 LITES 15-20 mmHg
- 0265 TruSheer 20-30 mmHg
- 0255 TruSheer 30-40 mmHg

Maternity Figure

- 1777 LITES 15-20 mmHg
- 0267 TruSheer 20-30 mmHg



1765 BLACK

Material Information

1765
Manufactured using 88% nylon and 12% spandex

1775 | 1777
Manufactured using 84% nylon and 16% spandex

0267
Manufactured using 71% nylon and 29% spandex

0265 | 0255
Manufactured using 57% nylon and 43% spandex


Product Description

Provide mild, moderate, firm, or extra firm graduated compression to help energize the legs, and reduce leg discomfort. They provide long-lasting support that appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm, or Extra-Firm support
- Fashionable sheer fabric
- Flat panty seams for a comfortable inconspicuous fit
- Ventilated cotton panel for comfort
- Available with maternity panel (1777 and 0267)

Available Colors



*Colors may not be available in all styles.
Please refer to the Product Summary on pages 38 and 39 for more information.


CIRCUMFERENCE

| | ANKLE | THIGH | HIPS | |
|---------|---------------|-----------------|-----------|----------------------|
| SMALL | 7" - 8.25" | 15.75" - 24.5" | 28" - 46" | 0265 0255 0267 |
| MEDIUM | 8.25" - 9.5" | 18" - 27.5" | 30" - 50" | |
| LARGE | 9.5" - 10.75" | 21.25" - 30.75" | 32" - 54" | |
| X-LARGE | 10.75" - 12" | 23.5" - 32" | 40" - 65" | |

Ladies' Sheer Compression Pantyhose

| Model | Compression | Description | Latex Free |
|-------|--------------|---------------------|------------|
| 1765 | 8 - 15 mmHg | Pantyhose | ✓ |
| 1775 | 15 - 20 mmHg | Pantyhose | ✓ |
| 1777 | 15 - 20 mmHg | Maternity Pantyhose | ✓ |
| 0265 | 20 - 30 mmHg | Pantyhose | ✓ |
| 0255 | 30 - 40 mmHg | Pantyhose | ✓ |
| 0267 | 20 - 30 mmHg | Maternity Pantyhose | ✓ |

WEIGHT



Ladies' Opaque

0373 15-20 mmHg
Knee High BLACK



Ladies' Opaque Knee High

Closed Toe

- 0373 15-20 mmHg
- 0363 20-30 mmHg

Open Toe

- 0371 15-20 mmHg
- 0361 20-30 mmHg



0373 BLACK

0361 BEIGE

Material Information

0371 | 0373 | 0361 | 0363
 Manufactured using 85% nylon and 15% spandex


Product Description

TRUFORM Opaque resemble high fashion hosiery yet feature moderate and firm graduated compression technology - very soft for the most comfortable stocking in its compression class. The "easy-on", "easy-off" construction is lightweight and highly effective in helping restore leg health. Long-lasting medical-grade support, fashionable look and super soft texture appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Moderate or Firm support
- Advanced microfiber-soft yarns and knit-in heel and toe offer comfort, better fit and durability
- Fashionably smooth and soft opaque fabric provides for year-round wear and covers blemishes
- Soft, stay-up top for all-day comfort

Available Colors



CIRCUMFERENCE

| | ANKLE | CALF |
|---------|---------------|----------------|
| SMALL | 7" - 8.25" | 11" - 15" |
| MEDIUM | 8.25" - 9.5" | 11.75" - 16.5" |
| LARGE | 9.5" - 10.75" | 12.5" - 18" |
| X-LARGE | 10.75" - 12" | 13" - 19" |

0373
0363
0371
0361

LATEX FREE

| Ladies' Opaque Knee High Compression Stockings | | | |
|--|--------------|--------------------|------------|
| Model | Compression | Description | Latex Free |
| 0371 | 15 - 20 mmHg | Knee High Open Toe | ✓ |
| 0373 | 15 - 20 mmHg | Knee High | ✓ |
| 0361 | 20 - 30 mmHg | Knee High Open Toe | ✓ |
| 0363 | 20 - 30 mmHg | Knee High | ✓ |

Ladies' Opaque Thigh High

Closed Toe

● 0364 20-30 mmHg

Open Toe

● 0362 20-30 mmHg

Silicone stay-up top



0364 BEIGE

0362 BLACK

Material Information

0364 | 0362
Manufactured using 85% nylon and 15% spandex

Product Description


TRUFORM Opaque resemble high fashion hosiery yet feature firm graduated compression technology - very soft for the most comfortable stocking in its compression class. The "easy-on", "easy-off" construction is lightweight and highly effective in helping restore leg health. Long-lasting medical-grade support, fashionable look and super soft texture appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Firm support
- Advanced microfiber-soft yarns and knit-in heel offer comfort, better fit and durability
- Fashionably smooth soft opaque fabric provides for year-round wear and helps conceal blemishes
- Functional silicone dot stay-up top for all-day wear

LATEX FREE

Available Colors



CIRCUMFERENCE

| | ANKLE | CALF | THIGH | |
|---------|---------------|----------------|-----------------|--------------|
| SMALL | 7" - 8.25" | 11" - 15" | 15.75" - 24.5" | 0364 0362 |
| MEDIUM | 8.25" - 9.5" | 11.75" - 16.5" | 18" - 27.5" | |
| LARGE | 9.5" - 10.75" | 12.5" - 18" | 21.25" - 30.75" | |
| X-LARGE | 10.75" - 12" | 13" - 19" | 23.5" - 32" | |

Ladies' Opaque Thigh High Compression Stockings

| Model | Compression | Description | Latex Free |
|-------|--------------|---------------------|------------|
| 0364 | 20 - 30 mmHg | Thigh High | ✓ |
| 0362 | 20 - 30 mmHg | Thigh High Open Toe | ✓ |

Ladies' Opaque Pantyhose

Closed Toe

● 0365 20-30 mmHg



0365 BLACK

0365 BIEGE

Material Information

365
Manufactured using 85% nylon and 15% spandex

Product Description


TRUFORM Opaque resemble high fashion hosiery yet feature firm graduated compression technology - very soft for the most comfortable stocking in its compression class. The "easy-on", "easy-off" construction is lightweight and highly effective in helping restore leg health. Long-lasting medical-grade support, fashionable look and super soft texture appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Firm support
- Advanced microfiber-soft yarns and knit-in heel and toe offer comfort, better fit and durability
- Fashionably smooth soft opaque fabric provides for year-round wear and helps conceal blemishes
- Flat panty seams for a comfortable inconspicuous fit
- Ventilated cotton panel for greater wearing comfort

LATEX FREE

Available Colors



CIRCUMFERENCE


| | ANKLE | THIGH | HIPS |
|---------|---------------|-----------------|-----------|
| SMALL | 7" - 8.25" | 15.75" - 24.5" | 28" - 46" |
| MEDIUM | 8.25" - 9.5" | 18" - 27.5" | 30" - 50" |
| LARGE | 9.5" - 10.75" | 21.25" - 30.75" | 32" - 54" |
| X-LARGE | 10.75" - 12" | 23.5" - 32" | 40" - 65" |

0365

Ladies' Opaque Compression Pantyhose

| Model | Compression | Description | Latex Free |
|-------|--------------|-------------|------------|
| 0365 | 20 - 30 mmHg | Pantyhose | ✓ |

Ladies' Pattern Socks



1976 15-20 mmHg
Knee High Diamond
BLACK




Ladies' Pattern Socks Knee High



1973 NAVY

1975 WHITE

1976 BLACK



Rib

- 1973 15-20 mmHg

Cable

- 1975 15-20 mmHg

Diamond

- 1976 15-20 mmHg

Material Information

1973 | 1975 | 1976

Manufactured using 85% nylon and 15% spandex

Product Description


TRUFORM trouser socks are designed to help energize the legs and reduce leg discomfort. Microfiber nylon on the rib, cable, and diamond patterns is soft, and cooler to wear than standard nylon. Attractive designer-knit delicate patterns make these socks an excellent choice for business, or business-casual wear.

Product Features

- Therapeutic graduated compression
- Moderate support
- Soft, comfortable microfiber nylon construction
- Smooth, seamless linked toe closure
- Extra-wide comfort band
- Attractive, designer-knit delicate patterns

LATEX FREE

Available Colors



CIRCUMFERENCE

| | ANKLE | CALF |
|---------|---------------|-------------|
| SMALL | 7" - 8.25" | 11" - 14" |
| MEDIUM | 8.25" - 9.75" | 13.5" - 16" |
| LARGE | 9.5" - 10.75" | 15.5" - 18" |
| X-LARGE | 11" - 12.5" | 17.5" - 20" |

1973
1975
1976

Ladies' Knee High Compression Socks

| Model | Compression | Description | Latex Free |
|-------|--------------|-------------------|------------|
| 1973 | 15 - 20 mmHg | Knee High Rib | ✓ |
| 1975 | 15 - 20 mmHg | Knee High Cable | ✓ |
| 1976 | 15 - 20 mmHg | Knee High Diamond | ✓ |

Men's Dress



1944 20-30 mmHg
Knee High BLACK

Men's Dress Knee | Thigh High

Knee High

- 1942 8-15 mmHg
- 1943 15-20 mmHg
- 1944 20-30 mmHg
- 1954 30-40 mmHg

Thigh High

- 1945 20-30 mmHg



Material Information

1942
Manufactured using 93% nylon and 7% spandex

1943 | 1944 | 1954 | 1945
Manufactured using 86% nylon and 14% spandex

Product Description


Truform Dress-style Support Socks for men provide mild, moderate, firm or extra firm graduated compression to help energize the legs - ideal when sitting or standing for long periods. Nylon-spandex materials are soft, and resemble specialty store, designer-knit socks. Contemporary rib-knit styling is suitable for business or business-casual settings.

Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm or Extra firm support available
- Roomy toe and defined heel for comfort and correct fit
- Contemporary rib-knit styling
- Functional, silicone dot stay-up top for all-day wear (1945)

LATEX FREE

Available Colors



*Colors may not be available in all styles.
Please refer to the Product Summary on pages 38 and 39 for more information.

Men's Knee/Thigh High Compression Socks

| Model | Compression | Description | Latex Free |
|-------|--------------|-------------|------------|
| 1942* | 8 - 15 mmHg | Knee High | ✓ |
| 1943 | 15 - 20 mmHg | Knee High | ✓ |
| 1944 | 20 - 30 mmHg | Knee High | ✓ |
| 1945 | 20 - 30 mmHg | Thigh High | ✓ |
| 1954 | 30 - 40 mmHg | Knee High | ✓ |

CIRCUMFERENCE

| | SHOE SIZE | ANKLE | CALF | MID-THIGH | LENGTH TO KNEE | |
|---------|-------------|---------------|-------------|-----------------|----------------|--------------------------------------|
| SMALL | 7 and under | 7" - 8.25" | 11" - 14" | 15.75" - 24.5" | up to 15" | 1942 1943 1944 1945 1954 |
| MEDIUM | 7 1/2 - 10 | 8.25" - 9.75" | 13.5" - 16" | 18" - 27.5" | up to 16" | |
| LARGE | 10 1/2 - 12 | 9.75" - 11" | 15.5" - 18" | 21.25" - 30.75" | up to 17" | |
| X-LARGE | 12 and over | 11" - 12.5" | 17.5" - 20" | 23.5" - 32" | up to 18" | |

*Use shoe size for 1942

Men's and Ladies' Casual



1933 15-20 mmHg
Casual Style Knee High BLACK



Men's & Ladies' Casual Socks

Ladies' Knee High

- 1963 15-20 mmHg

Men's Knee High

- 1933 15-20 mmHg
White color option available in athletic box

Men's Crew

- 1932 15-20 mmHg



1963 NAVY

1933 TAN

1932 WHITE

Material Information

1932 | 1933

Manufactured using 64% acrylic 29% nylon and 7% spandex

1963

Manufactured using 60% acrylic 33% nylon and 7% spandex

Product Description


Truform Men's & Ladies' Casual Support Socks are designed to help improve circulation and reduce leg discomfort. Shock-absorbing cushioned toe, sole and heel, moisture wicking fibers and anti-bacterial, stay-fresh finish provide all-day wearing comfort. Excellent for court sports, golfing, jogging, walking, or any leisure activity that requires comfort and support.

Product Features

- Therapeutic graduated compression
- Moderate support
- Shock-absorbing cushioned toe, sole and heel
- Moisture wicking acrylic materials
- Anti-bacterial, stay-fresh treatment helps control odor

LATEX FREE

Available Colors



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

Men's and Ladies' Casual Compression Socks

| Model | Compression | Description | Latex Free |
|----------------------|--------------|----------------------------|------------|
| 1963 Ladies' | 15 - 20 mmHg | Casual Style Knee High | ✓ |
| 1933 Men's (Colored) | 15 - 20 mmHg | Casual Style Knee High | ✓ |
| 1933 Men's (White) | 15 - 20 mmHg | Athletic Style Knee High | ✓ |
| 1932 Men's | 15 - 20 mmHg | Athletic Style Crew Length | ✓ |

| | CIRCUMFERENCE | | |
|---------|---------------|-------------|--------------|
| | ANKLE | CALF | |
| SMALL | 7" - 8.25" | 11" - 14" | 1933 1932 |
| MEDIUM | 8.25" - 9.75" | 13.5" - 16" | |
| LARGE | 9.75" - 11" | 15.5" - 18" | |
| X-LARGE | 11" - 12.5" | 17.5" - 20" | |

| | CIRCUMFERENCE | | |
|---------|---------------|----------------|------|
| | ANKLE | CALF | |
| SMALL | 7" - 8.25" | 11" - 15" | 1963 |
| MEDIUM | 8.25" - 9.5" | 11.75" - 16.5" | |
| LARGE | 9.5" - 10.75" | 12.5" - 18" | |
| X-LARGE | 10.75" - 12" | 13" - 19" | |

Medical Leg Wear



864 20-30 mmHg
Knee High Open Toe
BEIGE

Medical Leg Wear Below Knee

Below Knee Soft Top

- 8865 20-30 mmHg
- 8865S (short) 20-30 mmHg
- 8845 30-40 mmHg
- 8845S (short) 30-40 mmHg

Below Knee Dot Top

- 8864 20-30 mmHg
- 8844 30-40 mmHg

Silicone
dot top



8864 BEIGE



Material Information

8865(S) | 8845(S) | 8864 | 8844
Manufactured using 80% nylon and 20% spandex

Product Description

Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or Extra Firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top
- Standard and larger sizes

LATEX FREE

Available Colors



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

| | CIRCUMFERENCE | | LENGTH |
|----------|-------------------|---------------|---|
| | ANKLE | CALF | TO KNEE |
| SMALL | 7" - 8 1/4" | 11" - 14" | <p>Regular: Greater than 15" Short: 15" or less</p> |
| MEDIUM | 8 3/8" - 9 5/8" | 13 1/2" - 16" | |
| LARGE | 9 3/4" - 11" | 15 1/2" - 18" | |
| X-LARGE | 11 1/8" - 12 3/8" | 17 1/2" - 20" | |
| 2X-LARGE | 12 1/2" - 13 3/4" | 19 1/2" - 22" | |
| 3X-LARGE | 12 1/2" - 13 3/4" | 21 1/2" - 24" | up to 18" |

Classic Medical Knee High Compression Socks

| Model | Compression | Description | Latex Free |
|---------------|--------------|---------------------|------------|
| 8864 | 20 - 30 mmHg | Below Knee Dot Top | ✓ |
| 8865 | 20 - 30 mmHg | Below Knee Soft Top | ✓ |
| 8865S (Short) | 20 - 30 mmHg | Below Knee Soft Top | ✓ |
| 8844 | 30 - 40 mmHg | Below Knee Dot Top | ✓ |
| 8845 | 30 - 40 mmHg | Below Knee Soft Top | ✓ |
| 8845S (Short) | 30 - 40 mmHg | Below Knee Soft Top | ✓ |

Medical Leg Wear Below Knee Open Toe

Below Knee Open Toe Soft Top

- 0865 20-30 mmHg
- 0865S (short) 20-30 mmHg
- 0845 30-40 mmHg
- 0845S (short) 30-40 mmHg

Below Knee Open Toe Dot Top

- 0864 20-30 mmHg
- 0844 30-40 mmHg

Silicone
dot top



0864 BEIGE

Material Information

0865(S) | 0845(S) | 0864 | 0844
Manufactured using 80% nylon and 20% spandex

Product Description

Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. The open toe feature allows for easy toe inspection and permits one standard size to fit a variety of foot dimensions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or Extra Firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top
- Standard and larger sizes

LATEX FREE

Available Colors



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

| | CIRCUMFERENCE | | LENGTH |
|----------|-------------------|---------------|---|
| | ANKLE | CALF | TO KNEE |
| SMALL | 7" - 8 1/4" | 11" - 14" | Regular: Greater than 15" Short: 15" or less |
| MEDIUM | 8 3/8" - 9 5/8" | 13 1/2" - 16" | |
| LARGE | 9 3/4" - 11" | 15 1/2" - 18" | |
| X-LARGE | 11 1/8" - 12 3/8" | 17 1/2" - 20" | |
| 2X-LARGE | 12 1/2" - 13 3/4" | 19 1/2" - 22" | |
| 3X-LARGE | 12 1/2" - 13 3/4" | 21 1/2" - 24" | up to 18" |

Classic Medical Knee High Open Toe Compression Socks

| Model | Compression | Description | Latex Free |
|---------------|--------------|---------------------|------------|
| 0864 | 20 - 30 mmHg | Below Knee Dot Top | ✓ |
| 0865 | 20 - 30 mmHg | Below Knee Soft Top | ✓ |
| 0865S (Short) | 20 - 30 mmHg | Below Knee Soft Top | ✓ |
| 0844 | 30 - 40 mmHg | Below Knee Dot Top | ✓ |
| 0845 | 30 - 40 mmHg | Below Knee Soft Top | ✓ |
| 0845S (Short) | 30 - 40 mmHg | Below Knee Soft Top | ✓ |

Medical Leg Wear Thigh High

Thigh High Soft Top

- 8866 20-30 mmHg
- 8846 30-40 mmHg


Thigh High Dot Top

- 8868 20-30 mmHg
- 8848 30-40 mmHg

Thigh High Lace Top

- 8867 20-30 mmHg

Silicone lace
Stay-up top



Silicone dot
top



8867 BEIGE

8868 BEIGE



Material Information

8866 | 8868 | 8867
Manufactured using 80% nylon and 20% spandex

8846 | 8848
Manufactured using 83% nylon and 17% spandex

Product Description


Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or extra firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top

Note: A garter belt should be worn with soft top style numbers 8846 and 8866

Available Colors



*Colors may not be available in all styles.
Please refer to the Product Summary on pages 38 and 39 for more information.

LATEX FREE

| | CIRCUMFERENCE | | | LENGTH |
|---------|---|--|-----------|-----------|
| | ANKLE | CALF | THIGH | TO THIGH |
| SMALL | 7" - 8 ¹ / ₄ " | 11" - 14" | up to 20" | up to 26" |
| MEDIUM | 8 ³ / ₈ " - 9 ⁵ / ₈ " | 13 ¹ / ₂ " - 16" | up to 22" | up to 28" |
| LARGE | 9 ³ / ₄ " - 11" | 15 ¹ / ₂ " - 18" | up to 23" | up to 29" |
| X-LARGE | 11 ¹ / ₈ " - 12 ³ / ₈ " | 17 ¹ / ₂ " - 20" | up to 24" | up to 31" |

| Classic Medical Thigh High Compression Socks | | | |
|--|--------------|---------------------|------------|
| Model | Compression | Description | Latex Free |
| 8866 | 20 - 30 mmHg | Thigh High Soft Top | ✓ |
| 8867 | 20 - 30 mmHg | Thigh High Lace Top | ✓ |
| 8868 | 20 - 30 mmHg | Thigh High Dot Top | ✓ |
| 8846 | 30 - 40 mmHg | Thigh High Soft Top | ✓ |
| 8848 | 30 - 40 mmHg | Thigh High Dot Top | ✓ |

Medical Leg Wear Thigh High Open Toe

Thigh High Soft Top

- 0866 20-30 mmHg
- 0846 30-40 mmHg

Thigh High Dot Top

- 0868 20-30 mmHg
- 0848 30-40 mmHg

Silicone
dot top



0868 BEIGE

Material Information

0846
Manufactured using 80% nylon and 20% spandex

0866 | 0868 | 0848
Manufactured using 83% nylon and 17% spandex

Product Description

Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. The open toe feature allows for easy toe inspection and permits one standard size to fit a variety of foot dimensions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or extra firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top

Note: A garter belt should be worn with soft top style numbers 0846 and 0866

LATEX FREE

Available Colors

Beige
BG



| | CIRCUMFERENCE | | | LENGTH |
|---------|---|--|-----------|-----------|
| | ANKLE | CALF | THIGH | TO THIGH |
| SMALL | 7" - 8 ¹ / ₄ " | 11" - 14" | up to 20" | up to 26" |
| MEDIUM | 8 ³ / ₈ " - 9 ⁵ / ₈ " | 13 ¹ / ₂ " - 16" | up to 22" | up to 28" |
| LARGE | 9 ³ / ₄ " - 11" | 15 ¹ / ₂ " - 18" | up to 23" | up to 29" |
| X-LARGE | 11 ¹ / ₈ " - 12 ³ / ₈ " | 17 ¹ / ₂ " - 20" | up to 24" | up to 31" |

Classic Medical Thigh High Open Toe Compression Socks

| Model | Compression | Description | Latex Free |
|-------|--------------|---------------------|------------|
| 0866 | 20 - 30 mmHg | Thigh High Soft Top | ✓ |
| 0868 | 20 - 30 mmHg | Thigh High Dot Top | ✓ |
| 0846 | 30 - 40 mmHg | Thigh High Soft Top | ✓ |
| 0848 | 30 - 40 mmHg | Thigh High Dot Top | ✓ |

Medical Leg Wear Panty hose

Standard Figure

- 1756 20-30 mmHg

Maternity Figure

- 1757 20-30 mmHg



Full Figure

- 1758 20-30 mmHg



Material Information

1756 | 1757 | 1758
Manufactured using 80% nylon and 20% spandex

Product Description


Designed to help relieve moderate conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm support.
- Classic, concealing smooth-knit texture fabric
- Flat panty seams for a comfortable inconspicuous fit
- Ventilated cotton panel for comfort
- Available with maternity panel (1757)

LATEX FREE



Available Colors



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

| Classic Medical Compression Pantyhose | | | |
|---------------------------------------|--------------|------------------|------------|
| Model | Compression | Description | Latex Free |
| 1756 | 20 - 30 mmHg | Standard Figure | ✓ |
| 1757 | 20 - 30 mmHg | Maternity Figure | ✓ |
| 1758 | 20 - 30 mmHg | Full Figure | ✓ |



Anti-Embolism




Anti-Embolism Below Knee | Thigh Length

Below Knee

- 8808 | 8808S (Short) 18 mmHg

Below Knee Open Toe

- 0808 | 0808S (Short) 18 mmHg



Thigh Length

- 8810 18 mmHg

Thigh Length Open Toe

- 0810 18 mmHg



Material Information

8808 | 8808S | 0808 | 0808S
Manufactured using 89% nylon and 11% spandex

8810
Manufactured using 85% nylon and 15% spandex

0810
Manufactured using 82% nylon and 18% spandex

Product Description

Doctor-recommended for the recuperating patient, Truform anti-embolism Stockings are designed to help reduce the danger of thromboembolic disease for the reclining or bed-confined patient. They are also ideal for use with the post-surgical, home convalescing patient. The open toe feature allows for easy toe inspection by caregiver. Beige and black color options allow for a less clinical appearance during early ambulation.


Product Features

- Therapeutic graduated compression; 18 mmHg targeted at the ankle
- Non-constricting, soft stay-up top for comfort (8808, 0808)
- Open toe or closed toe designs
- Regular and Short lengths (Greater or less than 15" from heel to bend of knee)
- Pressure-relief panel in the upper thigh (8810, 0810)
- Silicone stay-up top (8810, 0810)

LATEX FREE

| | CIRCUMFERENCE | | | LENGTH | |
|----------|-------------------|---------------|-----------|---|-----------|
| | ANKLE | CALF | THIGH | TO KNEE | TO THIGH |
| SMALL | 7" - 8 1/4" | 11" - 14" | up to 21" | Regular: Greater than 15" Short: 15" or less | up to 28" |
| MEDIUM | 8 3/8" - 9 5/8" | 13 1/2" - 16" | up to 22" | | up to 29" |
| LARGE | 9 3/4" - 11" | 15 1/2" - 18" | up to 24" | | up to 30" |
| X-LARGE | 11 1/8" - 12 3/8" | 17 1/2" - 20" | up to 26" | | up to 32" |
| 2X-LARGE | 12 1/2" - 13 3/4" | 19 1/2" - 22" | N/A | up to 18" | N/A |
| 3X-LARGE | 12 1/2" - 13 3/4" | 21 1/2" - 24" | N/A | up to 18" | N/A |

Available Colors



*Colors may not be available in all styles.
Please refer to the Product Summary on pages 38 and 39 for more information.

| Anti-Embolism Knee and Thigh High Compression | | | |
|---|-------------|---------------------------|------------|
| Model | Compression | Description | Latex Free |
| 8808 | 18 mmHg | Below Knee | ✓ |
| 8808S (Short) | 18 mmHg | Below Knee Short | ✓ |
| 0808 | 18 mmHg | Below Knee Open Toe | ✓ |
| 0808S (Short) | 18 mmHg | Below Knee Open Toe Short | ✓ |
| 8810 | 18 mmHg | Thigh Length | ✓ |
| 0810 | 18 mmHg | Thigh Length Open Toe | ✓ |

Diabetic And Health Socks



16-20-30 mmHg
Knee High Open Toe
Grey



TruSoft Socks



- Mini-crew**
1911
- Crew**
● 1912 8-15 mmHg
- Knee High**
● 1913 8-15 mmHg

1913 BLACK

1912 BLACK

1911 WHITE

Material Information

1911 | 1912
Manufactured using 81% Acrylic 17% nylon and 2% spandex

1913
Manufactured using 75% Acrylic 22% nylon and 3% spandex

Product Description


TruSoft Socks surround the feet with full-cushioned comfort. Non-constricting, gentle, graduated compression helps keep sock in place to prevent bunching and sliding. Five gender-neutral sizes, moisture-wicking yarns, stay-fresh treatment, and NEW seamless linked toe closure help protect and comfort sensitive feet - ideal for the diabetic patient, or anyone with sensitive feet issues.

Product Features

- Ankle length, crew length, and knee high
- Mild graduated compression
- Moisture-wicking yarns
- Antibacterial, anti-fungal finish helps control odor
- NEW non-irritating seamless linked toe closure
- Full-cushioned foot

LATEX FREE

Available Colors



*Colors may not be available in all styles.
Please refer to the Product Summary on pages 38 and 39 for more information.

| | MALE SHOE SIZE | FEMALE SHOE SIZE | |
|---------|----------------|------------------|----------------------|
| X-SMALL | 3 1/2 - 5 | 4 1/2 - 6 1/2 | 1911 1912 1913 |
| SMALL | 6 - 8 | 7 - 9 | |
| MEDIUM | 8 1/2 - 10 | 9 1/2 - 11 | |
| LARGE | 10 1/2 - 12 | 11 1/2 - 13 | |
| X-LARGE | 12 1/2 - 14 | | |

| TruSoft Diabetic Socks | | | |
|------------------------|-------------|-------------|------------|
| Model | Compression | Description | Latex Free |
| 1911 | N/A | Mini Crew | ✓ |
| 1912 | 8 - 15 mmHg | Crew Length | ✓ |
| 1913 | 8 - 15 mmHg | Knee High | ✓ |

Comfort Care Diabetic Socks

Casual Crew
1914

Dress Crew
1915

Loose Fit Top Crew
1918



1918 BLACK

1915 BLACK

1914 WHITE

Material Information

1914 | 1915
Manufactured using 100% Cotton

Product Description


Comfort Care casual socks are designed for anyone with special needs in foot care. Non-constricting, gentle, elasticity helps keep sock in place to prevent bunching and sliding. The natural cotton materials helps reduce irritation and blistering and provides exceptional comfort, easy care and excellent absorbency. Suitable for men and women - ideal for the diabetic patient, or anyone with sensitive feet issues.

Product Features

- Crew length, cut extra-wide and extra-deep
- Seamless toe closure
- Non-binding and comfortable
- Can be worn over a cast
- Suitable for wear at night for warmth
- Excellent cover for use with medicinal or softening creams

LATEX FREE

Available Colors



| | SHOE SIZE |
|-----------|-----------|
| SMALL | 6 - 8½ |
| MED.-LGE. | 9 - 11 |
| X-LARGE | 10 - 14 |

1915
1914

| | MEN'S SHOE SIZE | LADIES' SHOE SIZE |
|--------|-----------------|-------------------|
| MEDIUM | 6 - 9 | 7 - 10 |
| LARGE | 9 - 12 | 10 - 13 |

1918

| Comfort Care Diabetic Socks | | | |
|-----------------------------|-------------|---------------------------|------------|
| Model | Compression | Description | Latex Free |
| 1914 | N/A | Casual Crew Length | ✓ |
| 1915 | N/A | Dress Crew Length | ✓ |
| 1918 | N/A | Loose Fit Top Crew Length | ✓ |

Upper Extremity

Ready-Wear Soft Top Sleeve

- 3315 15-20 mmHg
- 3325 20-30 mmHg

Ready-Wear Dot Top Sleeve

- 3316 15-20 mmHg
- 3326 20-30 mmHg

Ready-Wear Gauntlet

- 3312 15-20 mmHg
- 3322 20-30 mmHg



Material Information

3315 | 3325 | 3316 | 3326 | 3312 | 3322
 Manufactured using 73% nylon and 27% spandex

Product Description

Truform Ready To Wear Gauntlets and Sleeves are designed to help manage swelling following trauma, surgery, or radiation therapy. Lightweight microfiber knit construction allows for easy application, the soft air-permeable fabric provides extended wearing comfort. Available in three standard sizes and basic neutral colors.

Product Features

- Therapeutic, moderate and firm graduated compression
- Breathable, lightweight materials and construction is soft and easy to apply
- Silicone dot top to help sleeve remain in place
- For immediate wear with only 3 easy measurements
- Not made with natural rubber latex

Measuring:

For gauntlet

Measure around the palm and around the smallest part of the wrist as shown.


For arm sleeve

Measure around the smallest part of the wrist, around the mid-lower arm (largest part of the forearm), and around the mid-upper arm as shown.



LATEX FREE

Available Colors



*Colors may not be available in all styles.
 Please refer to the Product Summary on pages 38 and 39 for more information.

| Ready-Wear Compression Upper Extremity | | | |
|--|--------------|---------------------|------------|
| Model | Compression | Description | Latex Free |
| 3315 | 15 - 20 mmHg | Arm Sleeve Soft Top | ✓ |
| 3316 | 15 - 20 mmHg | Arm Sleeve Dot Top | ✓ |
| 3325 | 20 - 30 mmHg | Arm Sleeve Soft Top | ✓ |
| 3326 | 20 - 30 mmHg | Arm Sleeve Dot Top | ✓ |
| 3312 | 15 - 20 mmHg | Gauntlet | ✓ |
| 3322 | 20 - 30 mmHg | Gauntlet | ✓ |

CIRCUMFERENCE

| | PALM | WRIST | MID-LOWER | MID-UPPER | |
|--------|----------------|---------------|-----------------|-----------------|-------------------------------------|
| SMALL | 5.75" - 7.25" | 5.25" - 7.5" | 7" - 9" | 8.5" - 15.25" | 3315 3326 3325 3312 3316 3322 |
| MEDIUM | 7.5" - 8.5" | 6" - 8.5" | 9.25" - 11.25" | 10.5" - 17.25" | |
| LARGE | 8.75" - 10.25" | 6.75" - 9.25" | 11.25" - 13.25" | 12.75" - 19.75" | |



757 Donning
Gloves



752 Stocking Applicator

TRUFORM Accessories

757 Donning Gloves

The non-slip, easy grip surface on the palm and fingers help ease compression hosiery on and protect the fabric from being damaged by rings or fingernails.

750 Slip Sock

Often the most difficult part of applying medical compression stockings is getting them over the heel and ankle. The Slip Sock is designed to help open toe stockings slide smoothly onto the foot and over the ankle by reducing friction between the patient's skin and the compression stocking.

752 Stocking Applicator

The Truform stocking applicator allows donning of stockings with minimal bending and strain.



750 Slip Sock

TRUFORM Garter Belts

760 3" Garter Belt

762 2" Lightweight Garter Belt For Women

765 6" Garter Belt For Men

These standard style garter belts are practical and easy to wear. Garter attachments are completely adjustable in length and provide a convenient way to keep thigh high stockings in place. Suitable for both men and women.



760 3"
Garter Belt




762 2"
Garter Belt



765 6"
Garter Belt
for Men



HIGH PRECISION MANUFACTURING



TRUFORM[®]
LEG HEALTH

Knitting & Sewing Operations

Truform employs the latest methods of semi-automated manufacturing, and utilizes the best in technology to produce a world-class product at affordable prices - our prescribing physicians trust them, our many clients appreciate their extraordinary value.

Each machine is programmed to precise manufacturing specifications, and our technicians monitor every step of the knitting process. The combination of technology and knitting expertise ensures that Truform Stockings are among the highest quality, pressure gradient garments produced anywhere in the world.

We combine the finest materials and the most modern methods of 21st Century manufacturing, along with hand craftsmanship to produce Truform garments.

Quality Assurance

Quality assurance checks are performed throughout our entire manufacturing process to verify fabric integrity, fit and a therapeutic, pressure-gradient profile.



TRUFORM Product Summary

Anti-Embolism Stockings 18 mmHg pg. 31

| Model | Length | Toe Style | Color | SIZES |
|---------|--------------|------------|-------|-----------------------|
| 808BG | Below Knee | Open Toe | Beige | S-M-L-XL-2XL-3XL |
| 808WH | Below Knee | Open Toe | White | S-M-L-XL-2XL-3XL |
| 808SBL | Below Knee | Open Toe | Beige | Short Length S-M-L-XL |
| 810BG | Thigh Length | Open Toe | Beige | S-M-L-XL |
| 810WH | Thigh Length | Open Toe | White | S-M-L-XL |
| 8808BG | Below Knee | Closed Toe | Beige | S-M-L-XL-2XL-3XL |
| 8808BL | Below Knee | Closed Toe | Black | S-M-L-XL-2XL-3XL |
| 8808WH | Below Knee | Closed Toe | White | S-M-L-XL-2XL-3XL |
| 8808SBL | Below Knee | Closed Toe | Beige | Short Length S-M-L-XL |
| 8810BG | Thigh Length | Closed Toe | Beige | S-M-L-XL |
| 8810BL | Thigh Length | Closed Toe | Black | S-M-L-XL |
| 8810WH | Thigh Length | Closed Toe | White | S-M-L-XL |

Medical Stockings 20 - 30 mmHg pg. 25 -28

| Model | Length | Toe Style | Color | SIZES |
|---------|------------|------------|-------|--------------------------|
| 864BG | Below Knee | Open Toe | Beige | Dot Top S-M-L-XL-2XL-3XL |
| 865BG | Below Knee | Open Toe | Beige | S-M-L-XL-2XL-3XL |
| 865BL | Below Knee | Open Toe | Black | S-M-L-XL-2XL-3XL |
| 865WH | Below Knee | Open Toe | White | S-M-L-XL-2XL-3XL |
| 865SBL | Below Knee | Open Toe | Beige | Short Length S-M-L-XL |
| 865SBL | Below Knee | Open Toe | Black | Short Length S-M-L-XL |
| 866BG | Thigh High | Open Toe | Beige | S-M-L-XL |
| 868BG | Thigh High | Open Toe | Beige | Dot Top S-M-L-XL |
| 8864BG | Below Knee | Closed Toe | Beige | Dot Top S-M-L-XL-2XL-3XL |
| 8865BG | Below Knee | Closed Toe | Beige | S-M-L-XL-2XL-3XL |
| 8865SBL | Below Knee | Closed Toe | Black | Short Length S-M-L-XL |
| 8865BL | Below Knee | Closed Toe | Black | S-M-L-XL-2XL-3XL |
| 8865WH | Below Knee | Closed Toe | White | S-M-L-XL-2XL-3XL |
| 8866BG | Thigh High | Closed Toe | Beige | S-M-L-XL |
| 8867BG | Thigh High | Closed Toe | Beige | Lace Top S-M-L-XL |
| 8868BG | Thigh High | Closed Toe | Beige | Dot Top S-M-L-XL |
| 8868BL | Thigh High | Closed Toe | Black | Dot Top S-M-L-XL |

Medical Stockings 30 - 40 mmHg pg. 25 -28

| Model | Length | Toe Style | Color | SIZES |
|---------|------------|------------|-------|--------------------------|
| 844BG | Below Knee | Open Toe | Beige | Dot Top S-M-L-XL-2XL-3XL |
| 845BG | Below Knee | Open Toe | Beige | S-M-L-XL-2XL-3XL |
| 845BL | Below Knee | Open Toe | Black | S-M-L-XL-2XL-3XL |
| 845WH | Below Knee | Open Toe | White | S-M-L-XL-2XL-3XL |
| 845SBL | Below Knee | Open Toe | Beige | Short Length S-M-L-XL |
| 846BG | Thigh High | Open Toe | Beige | S-M-L-XL |
| 848BG | Thigh High | Open Toe | Beige | Dot Top S-M-L-XL |
| 8844BG | Below Knee | Closed Toe | Beige | Dot Top S-M-L-XL-2XL-3XL |
| 8845BG | Below Knee | Closed Toe | Beige | S-M-L-XL-2XL-3XL |
| 8845SBL | Below Knee | Closed Toe | Beige | Short Length S-M-L-XL |
| 8845BL | Below Knee | Closed Toe | Black | S-M-L-XL-2XL-3XL |
| 8845WH | Below Knee | Closed Toe | White | S-M-L-XL-2XL-3XL |
| 8846BG | Thigh High | Closed Toe | Beige | S-M-L-XL |
| 8848BG | Thigh High | Closed Toe | Beige | Dot Top S-M-L-XL |
| 8848BL | Thigh High | Closed Toe | Black | Dot Top S-M-L-XL |

Medical Pantyhose 20 - 30 mmHg pg. 29

| Model | Length | Toe Style | Color | SIZES |
|--------|-----------|------------|-------|----------------------|
| 1756BG | Pantyhose | Closed Toe | Beige | P-M-T-XT |
| 1756BL | Pantyhose | Closed Toe | Black | P-M-T-XT |
| 1757BG | Pantyhose | Closed Toe | Beige | Maternity P-M-T-XT |
| 1757BL | Pantyhose | Closed Toe | Black | Maternity P-M-T-XT |
| 1758BG | Pantyhose | Closed Toe | Beige | Full Figure P-M-T-XT |
| 1758BL | Pantyhose | Closed Toe | Black | Full Figure P-M-T-XT |

Truform Lites Sheer Support Hosiery pg. 9 - 13

| Model | Length | Compression | Color | SIZES |
|--------|------------|--------------|-------------|--------------------------|
| 1763BG | Knee High | 8 - 15 mmHg | Beige | S-M-L-XL |
| 1763BL | Knee High | 8 - 15 mmHg | Black | S-M-L-XL |
| 1763NV | Knee High | 8 - 15 mmHg | Navy | S-M-L-XL |
| 1763TP | Knee High | 8 - 15 mmHg | Taupe | S-M-L-XL |
| 1764BG | Thigh High | 8 - 15 mmHg | Beige | S-M-L-XL |
| 1764BL | Thigh High | 8 - 15 mmHg | Black | S-M-L-XL |
| 1765BG | Pantyhose | 8 - 15 mmHg | Beige | P-M-T-XT-Q-QP |
| 1765BL | Pantyhose | 8 - 15 mmHg | Black | P-M-T-XT-Q-QP |
| 1765TP | Pantyhose | 8 - 15 mmHg | Taupe | P-M-T-XT-Q-QP |
| 1772ND | Knee High | 15 - 20 mmHg | Nude | Open Toe S-M-L-XL |
| 1772BL | Knee High | 15 - 20 mmHg | Black | Open Toe S-M-L-XL |
| 1773BG | Knee High | 15 - 20 mmHg | Beige | S-M-L-XL-2XL-3XL |
| 1773BL | Knee High | 15 - 20 mmHg | Black | S-M-L-XL |
| 1773IV | Knee High | 15 - 20 mmHg | Ivory | S-M-L-XL |
| 1773LB | Knee High | 15 - 20 mmHg | Light Beige | S-M-L-XL |
| 1773ND | Knee High | 15 - 20 mmHg | Nude | S-M-L-XL |
| 1773NV | Knee High | 15 - 20 mmHg | Navy | S-M-L-XL |
| 1773TP | Knee High | 15 - 20 mmHg | Taupe | S-M-L-XL |
| 1773WH | Knee High | 15 - 20 mmHg | White | S-M-L-XL |
| 1774BG | Thigh High | 15 - 20 mmHg | Beige | S-M-L-XL |
| 1774BL | Thigh High | 15 - 20 mmHg | Black | S-M-L-XL |
| 1774IV | Thigh High | 15 - 20 mmHg | Ivory | S-M-L-XL |
| 1774ND | Thigh High | 15 - 20 mmHg | Nude | S-M-L-XL |
| 1774TP | Thigh High | 15 - 20 mmHg | Taupe | S-M-L-XL |
| 1774WH | Thigh High | 15 - 20 mmHg | White | S-M-L-XL |
| 1775BG | Pantyhose | 15 - 20 mmHg | Beige | P-M-T-XT-Q-QP |
| 1775BL | Pantyhose | 15 - 20 mmHg | Black | P-M-T-XT-Q-QP |
| 1775IV | Pantyhose | 15 - 20 mmHg | Ivory | P-M-T-XT-Q-QP |
| 1775LB | Pantyhose | 15 - 20 mmHg | Light Beige | P-M-T-XT-Q-QP |
| 1775ND | Pantyhose | 15 - 20 mmHg | Nude | P-M-T-XT-Q-QP |
| 1775NV | Pantyhose | 15 - 20 mmHg | Navy | P-M-T-XT-Q-QP |
| 1775TP | Pantyhose | 15 - 20 mmHg | Taupe | P-M-T-XT-Q-QP |
| 1775WH | Pantyhose | 15 - 20 mmHg | White | P-M-T-XT-Q-QP |
| 1777BG | Pantyhose | 15 - 20 mmHg | Beige | Maternity P-M-T-XT |
| 1777BL | Pantyhose | 15 - 20 mmHg | Black | Maternity P-M-T-XT |
| 1782BL | Knee High | 15 - 20 mmHg | Black | Dot Pattern S-M-L-XL |
| 1782CH | Knee High | 15 - 20 mmHg | Charcoal | Dot Pattern S-M-L-XL |
| 1782ES | Knee High | 15 - 20 mmHg | Espresso | Dot Pattern S-M-L-XL |
| 1783BL | Knee High | 15 - 20 mmHg | Black | Diamond Pattern S-M-L-XL |
| 1783CH | Knee High | 15 - 20 mmHg | Charcoal | Diamond Pattern S-M-L-XL |
| 1783ES | Knee High | 15 - 20 mmHg | Espresso | Diamond Pattern S-M-L-XL |

Truform TruSheer Therapeutic Support Hosiery pg. 9 - 13

| Model | Length | Compression | Color | SIZES |
|-------|------------|--------------|-------|--------------------|
| 253BG | Knee High | 30 - 40 mmHg | Beige | S-M-L-XL |
| 253BL | Knee High | 30 - 40 mmHg | Black | S-M-L-XL |
| 253TP | Knee High | 30 - 40 mmHg | Taupe | S-M-L-XL |
| 254BG | Thigh High | 30 - 40 mmHg | Beige | S-M-L-XL |
| 254BL | Thigh High | 30 - 40 mmHg | Black | S-M-L-XL |
| 254TP | Thigh High | 30 - 40 mmHg | Taupe | S-M-L-XL |
| 255BG | Pantyhose | 30 - 40 mmHg | Beige | S-M-L-XL |
| 255BL | Pantyhose | 30 - 40 mmHg | Black | S-M-L-XL |
| 255TP | Pantyhose | 30 - 40 mmHg | Taupe | S-M-L-XL |
| 263BG | Knee High | 20 - 30 mmHg | Beige | S-M-L-XL |
| 263BL | Knee High | 20 - 30 mmHg | Black | S-M-L-XL |
| 263TP | Knee High | 20 - 30 mmHg | Taupe | S-M-L-XL |
| 264BG | Thigh High | 20 - 30 mmHg | Beige | S-M-L-XL |
| 264BL | Thigh High | 20 - 30 mmHg | Black | S-M-L-XL |
| 264TP | Thigh High | 20 - 30 mmHg | Taupe | S-M-L-XL |
| 265BG | Pantyhose | 20 - 30 mmHg | Beige | S-M-L-XL |
| 265BL | Pantyhose | 20 - 30 mmHg | Black | S-M-L-XL |
| 265TP | Pantyhose | 20 - 30 mmHg | Taupe | S-M-L-XL |
| 267BG | Pantyhose | 20 - 30 mmHg | Beige | Maternity S-M-L-XL |
| 267BL | Pantyhose | 20 - 30 mmHg | Black | Maternity S-M-L-XL |

Truform Ladies' Therapeutic Support Socks pg. 19 & 23

| Model | Length | Compression | Color | SIZES |
|--------|-----------|--------------|-------|--------------------------|
| 1963BL | Knee High | 15 - 20 mmHg | Black | Casual Style S-M-L-XL |
| 1963NV | Knee High | 15 - 20 mmHg | Navy | Casual Style S-M-L-XL |
| 1963TN | Knee High | 15 - 20 mmHg | Tan | Casual Style S-M-L-XL |
| 1963WH | Knee High | 15 - 20 mmHg | White | Casual Style S-M-L-XL |
| 1973BL | Knee High | 15 - 20 mmHg | Black | Rib Pattern S-M-L-XL |
| 1973BN | Knee High | 15 - 20 mmHg | Brown | Rib Pattern S-M-L-XL |
| 1973NV | Knee High | 15 - 20 mmHg | Navy | Rib Pattern S-M-L-XL |
| 1973TN | Knee High | 15 - 20 mmHg | Tan | Rib Pattern S-M-L-XL |
| 1973WH | Knee High | 15 - 20 mmHg | White | Rib Pattern S-M-L-XL |
| 1975BL | Knee High | 15 - 20 mmHg | Black | Cable Pattern S-M-L-XL |
| 1975BN | Knee High | 15 - 20 mmHg | Brown | Cable Pattern S-M-L-XL |
| 1975NV | Knee High | 15 - 20 mmHg | Navy | Cable Pattern S-M-L-XL |
| 1975TN | Knee High | 15 - 20 mmHg | Tan | Cable Pattern S-M-L-XL |
| 1975WH | Knee High | 15 - 20 mmHg | White | Cable Pattern S-M-L-XL |
| 1976BL | Knee High | 15 - 20 mmHg | Black | Diamond Pattern S-M-L-XL |
| 1976BN | Knee High | 15 - 20 mmHg | Brown | Diamond Pattern S-M-L-XL |
| 1976NV | Knee High | 15 - 20 mmHg | Navy | Diamond Pattern S-M-L-XL |
| 1976TN | Knee High | 15 - 20 mmHg | Tan | Diamond Pattern S-M-L-XL |
| 1976WH | Knee High | 15 - 20 mmHg | White | Diamond Pattern S-M-L-XL |

Ladies' Opaque Therapeutic Support Hosiery pg. 15 - 17

| Model | Length | Compression | Color | SIZES |
|--------|------------|--------------|-------|-------------------|
| 0371BG | Knee High | 15 - 20 mmHg | Beige | Open Toe S-M-L-XL |
| 0371BL | Knee High | 15 - 20 mmHg | Black | Open Toe S-M-L-XL |
| 0373BG | Knee High | 15 - 20 mmHg | Beige | S-M-L-XL |
| 0373BL | Knee High | 15 - 20 mmHg | Black | S-M-L-XL |
| 0361BG | Knee High | 20 - 30 mmHg | Beige | Open Toe S-M-L-XL |
| 0361BL | Knee High | 20 - 30 mmHg | Black | Open Toe S-M-L-XL |
| 0363BG | Knee High | 20 - 30 mmHg | Beige | S-M-L-XL |
| 0363BL | Knee High | 20 - 30 mmHg | Black | S-M-L-XL |
| 0362BG | Knee High | 20 - 30 mmHg | Beige | Open Toe S-M-L-XL |
| 0362BL | Thigh High | 20 - 30 mmHg | Black | Open Toe S-M-L-XL |
| 0364BG | Thigh High | 20 - 30 mmHg | Beige | S-M-L-XL |
| 0364BL | Thigh High | 20 - 30 mmHg | Black | S-M-L-XL |
| 0365BG | Pantyhose | 20 - 30 mmHg | Beige | S-M-L-XL |
| 0365BL | Pantyhose | 20 - 30 mmHg | Black | S-M-L-XL |

Truform Men's Therapeutic Support Socks pg. 21 & 23

| Model | Length | Compression | Color | SIZES |
|--------|-------------|--------------|-------|-------------------------|
| 1932WH | Crew Length | 15 - 20 mmHg | White | Athletic Style S-M-L-XL |
| 1933BL | Knee High | 15 - 20 mmHg | Black | Casual Style S-M-L-XL |
| 1933BN | Knee High | 15 - 20 mmHg | Brown | Casual Style S-M-L-XL |
| 1933NV | Knee High | 15 - 20 mmHg | Navy | Casual Style S-M-L-XL |
| 1933TN | Knee High | 15 - 20 mmHg | Tan | Casual Style S-M-L-XL |
| 1933WH | Knee High | 15 - 20 mmHg | White | Athletic Style S-M-L-XL |
| 1942BL | Knee High | 8 - 15 mmHg | Black | Dress Style S-M-L-XL |
| 1942BN | Knee High | 8 - 15 mmHg | Brown | Dress Style S-M-L-XL |
| 1942NV | Knee High | 8 - 15 mmHg | Navy | Dress Style S-M-L-XL |
| 1942TN | Knee High | 8 - 15 mmHg | Tan | Dress Style S-M-L-XL |
| 1943BL | Knee High | 15 - 20 mmHg | Black | Dress Style S-M-L-XL |
| 1943BN | Knee High | 15 - 20 mmHg | Brown | Dress Style S-M-L-XL |
| 1943NV | Knee High | 15 - 20 mmHg | Navy | Dress Style S-M-L-XL |
| 1943TN | Knee High | 15 - 20 mmHg | Tan | Dress Style S-M-L-XL |
| 1943WH | Knee High | 15 - 20 mmHg | White | Dress Style S-M-L-XL |
| 1944BL | Knee High | 20 - 30 mmHg | Black | Dress Style S-M-L-XL |
| 1944BN | Knee High | 20 - 30 mmHg | Brown | Dress Style S-M-L-XL |
| 1944NV | Knee High | 20 - 30 mmHg | Navy | Dress Style S-M-L-XL |
| 1944TN | Knee High | 20 - 30 mmHg | Tan | Dress Style S-M-L-XL |
| 1945BL | Thigh High | 20 - 30 mmHg | Black | Dress Style S-M-L-XL |
| 1954BL | Knee High | 30 - 40 mmHg | Black | Dress Style S-M-L-XL |
| 1954BN | Knee High | 30 - 40 mmHg | Brown | Dress Style S-M-L-XL |
| 1954NV | Knee High | 30 - 40 mmHg | Navy | Dress Style S-M-L-XL |
| 1954TN | Knee High | 30 - 40 mmHg | Tan | Dress Style S-M-L-XL |

Health & Diabetic TruSoft & Comfort Care Socks pg. 33 - 34

| Model | Length | Compression | Color | SIZES |
|--------|-------------|-------------|-------|----------------------|
| 1911WH | Mini-Crew | N/A | White | TruSoft XS-S-M-L-XL |
| 1912BL | Crew Length | 8 - 15 mmHg | Black | TruSoft XS-S-M-L-XL |
| 1912WH | Crew Length | 8 - 15 mmHg | White | TruSoft XS-S-M-L-XL |
| 1913BL | Knee High | 8 - 15 mmHg | Black | TruSoft XS-S-M-L-XL |
| 1913WH | Knee High | 8 - 15 mmHg | White | TruSoft XS-S-M-L-XL |
| 1914WH | Crew Length | N/A | White | Casual Socks S-ML-XL |
| 1914BL | Crew Length | N/A | Black | Casual Socks S-ML-XL |
| 1915WH | Crew Length | N/A | White | Dress Socks S-ML-XL |
| 1915BL | Crew Length | N/A | Black | Dress Socks S-ML-XL |
| 1918BL | Crew Length | N/A | Black | Loose Fit Top M-L |
| 1918WH | Crew Length | N/A | White | Loose Fit Top M-L |

COMPRESSION
THERAPY COLLECTION



TRUFORM is a Division of SAI Therapeutic Brands
ph: 800-888-0458 • fax: 800-309-9055 • truform.com

3960 Rosslyn Dr, Cincinnati, OH
45209, United States

In Canada, a Division of Airway Surgical Appliances Ltd.
ph: 800-267-3476 • fax: 800-267-7805 • truform.com

189 Colonnade Rd S, Ottawa, ON
K2E 7J4, Canada



TRUFORM[®]
LEG HEALTH